*Instructions for use
Istruzioni per l'uso
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Указания за употреба
Upute za upotrebu
Пайдалану бойынша нұсқаулар
Instruções de utilização*

*Instrucțiuni de utilizare
Käyttöohjeet*

R9900101 | FLOW REVERSE

CE
0051

EN	<i>Accessory for Hyperthermic Perfusion</i>
IT	<i>Accessorio per Perfusione Ipertermica</i>
FR	<i>Accessoire pour Perfusion Hyperthermique</i>
DE	<i>Zubehör für hypertherme Perfusion</i>
ES	<i>Accesorio para Perfusión Hipertérmica</i>
BG	<i>Акcesoар за хипертермична перфузия</i>
HR	<i>Dodatak za hipertermijsku perfuziju</i>
KK	<i>Гипертермикалық перфузия үшін қосалқы құрал</i>
PT	<i>Acessório para Perfusão Hipertérmica</i>
RO	<i>Accesoriu pentru perfuzie hipertermică</i>
FI	<i>Lisävarusteet hypertermiseen perfuusioon</i>

EN. Do not reuse
IT. Non riutilizzare
FR. Ne pas réutiliser
DE. Nicht zur Wiederverwendung
ES. No reutilizar
TR. Tek kullanımlık
NL. Niet opnieuw gebruiken
RU. Не использовать повторно
PL. Do jednorazowego użytku
SR. Nije za višekratnu upotrebu
LT. Nenaudokite pakartotinai
SV. Får inte återanvändas
BG. Да не се използва повторно
DA. Må ikke genbruges
HR. Nije za višekratnu upotrebu
ZH. 请不要重复使用
ZH. 不可重複使用
CS. Nepoužívejte opakovaně
KK. Қайта пайдалануға болмайды
PT. Não reutilizar
RO. A nu se reutiliza
AR. لا تُعد استخداماً
FI. Ei saa käyttää uudelleen

EN. Content
IT. Contenuto
FR. Contenu
DE. Inhalt
ES. Contenido
TR. İçerik
NL. Inhoud
RU. Содержит
PL. Zawiera
SR. Sadržaj
LT. Sudėtis
SV. Innehåll
BG. Съдържание
DA. Indhold
HR. Sadržaj
ZH. 内容
ZH. 潤德
CS. Obsah
KK. Содержимое
PT. Conteúdo
RO. Conținut
AR. المحتوى
FI. Sisältö

EN. Non-pyrogenic
IT. Apirogeno
FR. Apyrogène
DE. Pyrogenfrei
ES. Apirogeno
TR. Apirojen
NL. Niet-pyrogeen
RU. Апирогенно
PL. Niepirogenny
SR. Ne stvara toplotu
LT. Nepirogeniškas
SV. Pyrogenfri
BG. Апирогенно
DA. Ikke-pyrogen
HR. Nepirogeno
ZH. 无热源的
ZH. 無致熱原
CS. Nepyrogní
KK. Апирогенді
PT. Não pirogénico
RO. Apirogen
AR. غير بروجيني
FI. Ei-pyrogeeninen

EN. Quantity
IT. Quantità
FR. Quantité
DE. Menge
ES. Cantidad
TR. Adet
NL. Aantal
RU. Количество
PL. Ilość
SR. Količina
LT. Kiekis
SV. Antal
BG. Количество
DA. Antal
HR. Količina
ZH. 数量
ZH. 數量
CS. Množství
KK. Саны
PT. Quantidade
RO. Cantitate

EN. Keep dry
IT. Teme l'umidità
FR. Craint l'humidité
DE. Trocken lagern
ES. Matener seco
TR. Kuru saklayın
NL. Droog bewaren
RU. Держать сухим
PL. Przechowywać w suchym miejscu
SR. Čuvati na suvom
LT. Laikykite sausai
SV. Förvaras torrt
BG. Да се държи на сухо място
DA. Skal holdes tørt
HR. Čuvati na suhom
ZH. 请保持干燥
ZH. 保持乾燥
CS. Uchovávejte v suchu
KK. Құрғақ сақтаңыз
PT. Manter seco
RO. A se păstra ferit de umezeală
AR. حافظ على جفافه
FI. Säilytettävä kuivana

EN. Date of manufacture
IT. Data di fabbricazione
FR. Date de fabrication
DE. Herstellungsdatum
ES. Fecha de fabricación
TR. Üretim tarihi
NL. Fabricagedatum
RU. Дата производства
PL. Data produkcji
SR. Datum proizvodnje
LT. Pagaminimo data
SV. Tillverkningsdatum
BG. Дата на производство
DA. Fremstillingsdato
HR. Datum proizvodnje
ZH. 生产日期
ZH. 製造日期
CS. Datum výroby
KK. Өндіру күні
PT. Data de fabrico
RO. Data fabricației
AR. تاريخ التصنيع
FI. Valmistuspäivämäärä

EN. Do not re-sterilize
IT. Non ristilizzare
FR. Ne pas restérilizer
DE. Nicht reesterilisieren
ES. No volver a esterilizar
TR. Tekrar steril etmeyin
NL. Niet opnieuw steriliseren
RU. Не стерилизовать повторно
PL. Nie sterylizować
SR. Ne sterilizovati više puta
LT. Nesterilizuoti pakartotinai
SV. Får inte återsteriliseras
BG. Да не се стерилизира повторно
DA. Må ikke gensteriliseres
HR. Ponovna sterilizacija nije dopuštena
ZH. 请勿重复消毒
ZH. 不可重複滅菌
CS. Nesterilizujte opakovaně
KK. Қайта стерилдеуге болмайды
PT. Não reesterilizar
RO. A nu se reesteriliza
AR. لا تُعد تعقيم
FI. Ei saa steriloida uudelleen

EN. Use by
IT. Utilizzare entro
FR. Utiliser jusqu'à
DE. Verwendbar bis
ES. Fecha de caducidad
TR. Son kullanım tarihi
NL. Te gebruiken vóór
RU. Использовать до
PL. Zużyć do
SR. Rok upotrebe
LT. Tinka naudoti iki
SV. Använd före
BG. Годен до
DA. Anvendes inden
HR. Upotrijebiti do
ZH. 使用者
ZH. 保存期限
CS. Spotřebujte do
KK. Жарамдылық мерзімі
PT. Data de validade
RO. A se utiliza înainte de

EN. Batch code
IT. Numero di lotto
FR. Numéro de lot
DE. Chargenbezeichnung
ES. Número de lote
TR. Lot numarası
NL. Partijnummer
RU. Номер партии
PL. Kod partii
SR. Šifra serije
LT. Partijos numeris
SV. Batchnummer
BG. Партиден номер
DA. Partikode
HR. Šifra serije
ZH. 生产批号
ZH. 製造批號
CS. Kód šarže
KK. Жинақ коды
PT. Código de lote
RO. Număr lot
AR. رمز التشغيل
FI. Eräkkoodi

EN. Manufacturer
IT. Fabbricante
FR. Fabricant
DE. Hersteller
ES. Fabricante
TR. Üretici
NL. Fabrikant
RU. Производитель
PL. Producent
SR. Proizvođač
LT. Gamintojas
SV. Tillverkare
BG. Производител
DA. Producent
HR. Proizvođač
ZH. 生产厂家
ZH. 製造廠
CS. Výrobce
KK. Өндіруші
PT. Fabricante
RO. Producător
AR. الشركة المصنعة
FI. Valmistaja

EN. Assembled by
IT. Assemblato da
FR. Assemblé par
DE. Zusammengefügt von
ES. Ensamblado por
TR. Birleştirici
NL. Geassembleerd door
RU. Упаковщик
PL. Złożono
SR. Sastavio
LT. Surinko
SV. Monterad av
BG. Сглобено от
DA. Samlet af
HR. Sastavio
ZH. 组装
CS. Sestavil
KK. Құрастырушы
PT. Montado por
RO. Asamblat de
AR. التجميع من قبل
FI. Kokoonpanija

REF

AR.	الكمية
FI.	Määrä
EN.	Catalogue number
IT.	Numero di catalogo
FR.	Référence du catalogue
DE.	Bestellnummer
ES.	Número de catálogo
TR.	Referans numarası
NL.	Catalogusnummer
RU.	Номер по каталогу
PL.	Numer katalogowy
SR.	Kataloški broj
LT.	Katalogo numeris
SV.	Katalognummer
BG.	Каталожен номер
DA.	Katalognummer
HR.	Broj u katalogu
ZH.	目录号
ZH.	產品型號
CS.	Katalogové číslo
KK.	Каталог нөмірі
PT.	Número de catálogo
RO.	Număr de catalog
AR.	رقم الكتالوج
FI.	Luettelonumero

EN.	Keep away from heat
IT.	Tenere lontano da sorgenti di calore
FR.	Craint la chaleur
DE.	Vor Wärmequellen fernhalten
ES.	Mantener lejos de fuentes de calor
TR.	Isı kaynaklarından uzak tutun
NL.	Uit de buurt van warmtebronnen houden
RU.	Беречь от тепла
PL.	Przechowywać z dala od źródeł ciepła
SR.	Držati dalje od toplote
LT.	Saugokite nuo šilumos
SV.	Skydda från värme
BG.	Да се пази от топлина
DA.	Må ikke udsættes for varme
HR.	Držati dalje od topline
ZH.	远离高温
ZH.	避免日照
CS.	Chraňte před teplem
KK.	Ыстыктан алыс ұстаңыз
PT.	Manter afastado do calor
RO.	A se păstra departe de sursele de căldură
AR.	ابتقه بعيداً عن الحرارة
FI.	Suojattava lämmöltä

STERILE EO

EN.	Sterile - Ethylene oxide sterilized
IT.	Sterile - Sterilizzato ad ossido d'etilene
FR.	Stérile - Stérilisé à l'oxyde d'éthylène
DE.	Steril - Mit Äthylenoxid sterilisiert
ES.	Estéril - Esterilizado por óxido de etileno
TR.	Steril - Etilen Oksit
NL.	Steriel - Gesteriliseerd met ethyleenoxide
RU.	Стерильно - Стерилизовано окисью этилена
PL.	Sterylnie - Wysterylizowane tlenkiem etylenu
SR.	Sterilno - Sterilisano etilen oksidom
LT.	Sterilu - Sterilizuota etileno oksidu
SV.	Steril - Steriliserad med etylenoxid
BG.	Стерилно - Стерилизирано с этиленов оксид
DA.	Steril - Steriliseret med ethylenoxid
HR.	Sterilno - Sterilizirano etilen-oksdom
ZH.	消毒 - 环氧乙烷消毒法
ZH.	EO滅菌
CS.	Sterilní - Sterilizováno ethylenoxidem
KK.	Стерильді - Этилен тотығымен стерильденген
PT.	Estéril - Esterilizado por óxido de etileno
RO.	Steril - Sterilizat cu oxid de etilenă
AR.	معقم - معقم بأكسيد الإيثيلين
FI.	Steriliili - steriloitu etyleenioksidilla

EN.	Not made with natural rubber latex
IT.	Non fatto con lattice di gomma naturale
FR.	Latex exempt
DE.	Latex-frei
ES.	Libre de látex
TR.	Lateks içermez
NL.	Latexvrij
RU.	Не содержит латекса
PL.	Nie zawiera lateksu
SR.	Ne sadrži lateks
LT.	Be latekso
SV.	Latexfri
BG.	Не съдържа латекс
DA.	Latex-fri
HR.	Bez lateksa
ZH.	请与乳胶隔离
ZH.	非乳膠製品
CS.	Neobsahuje latex
KK.	Латексціз

AR.	يستخدم من قبل
FI.	Viimeinen käyttöpäivä
EN.	Temperature limitation
IT.	Limiti di temperatura
FR.	Limites de température
DE.	Zulässiger Temperaturbereich
ES.	Límite de temperatura
TR.	Sıcaklık limitleri
NL.	Temperatuurlimieten
RU.	Пределы температур
PL.	Ograniczenia temperatury
SR.	Temperaturno ograničenje
LT.	Temperatūros apribojimas
SV.	Temperaturbegränsning
BG.	Температурно ограничение
DA.	Temperaturbegrænsning
HR.	Ograničenje temperature
ZH.	温度限制
ZH.	溫度限制
CS.	Teplotní omezení
KK.	Температура шери
PT.	Limite de temperatura
RO.	Limita de temperatură
AR.	حدود درجة الحرارة
FI.	Lämpötilarajoitus

EN.	Caution, consult accompanying documents
IT.	Attenzione, vedere le istruzioni per l'uso
FR.	Attention voir notice d'instructions
DE.	Achtung, Begleitdokumente beachten
ES.	Atención, ver instrucciones de uso
TR.	Dikkat, kullanım talimatlarına bakın
NL.	Attentie, de bijgevoegde documenten raadplegen
RU.	Внимание! Ознакомиться с сопроводительными документами
PL.	Uwaga! Zapoznać się z dołączonymi dokumentami
SR.	Pažnja, pogledajte prateću dokumentaciju
LT.	Dėmesio, perskaitykite pridedamus dokumentus
SV.	Varning – läs de medföljande dokumenten
BG.	Внимание, напратете справка в придружаващата документация
DA.	Forsigtig, rådfør dig med de medfølgende dokumenter
HR.	Oprez, pogledati prateće dokumente
ZH.	注意. 请阅读所附文件
ZH.	注意事項，請參考隨附文件
CS.	Upozornění: prostudujte si průvodní dokumentaci
KK.	Ескерту, бірге келген құжаттарды қараңыз
PT.	Cuidado, consultar os documentos fornecidos
RO.	Atenție, consultați documentele însoțitoare
AR.	احترس، راجع الوثائق المرفقة
FI.	Huomio, katso liiteasiakirjat

EN.	Do not use if package is damaged
IT.	Non utilizzare se la confezione è danneggiata
FR.	Ne pas utiliser si l'emballage est endommagé
DE.	Inhalt beschädigter Packung nicht verwenden
ES.	No usar si el paquete está dañado
TR.	Paket hasar görmüşse kullanmayın
NL.	Niet gebruiken als de verpakking beschadigd is
RU.	Не использовать, если упаковка повреждена
PL.	Nie używać, jeśli opakowanie jest uszkodzone
SR.	Ne koristiti ukoliko je pakovanje oštećeno
LT.	Nenaudokite, jeigu pakuota yra pažeista
SV.	Använd inte om förpackningen är skadad
BG.	Да не се използва, ако опаковката е повредена
DA.	Må ikke anvendes, hvis emballagen er beskadiget
HR.	Ne upotrebljavati ako je pakiranje oštećeno
ZH.	包装出现损坏的情况下请不要使用
ZH.	包裝破損不可使用
CS.	Nepoužívejte, pokud je poškozený obal
KK.	Орам зақымдалған болса, пайдаланбаңыз
PT.	Não utilizar se a embalagem estiver danificada
RO.	A nu se utiliza dacă ambalajul este deteriorat
AR.	لا تستخدمه في حالة تلف الحزمة
FI.	Älä käyttää, jos pakkaus on vaurioitunut

MD

EN.	Medical device
IT.	Dispositivo medico
FR.	Dispositif médical
DE.	Medizinprodukt
ES.	Producto sanitario
TR.	Tıbbi cihaz
NL.	Medisch hulpmiddel
RU.	Медицинское изделие
PL.	Wyrób medyczny
SR.	Medicinsko sredstvo
LT.	Medicinos priemonė
SV.	Medicinteknisk produkt
BG.	Медицинско изделие
DA.	Medicinsk udstyr
HR.	Medicinski proizvod
ZH.	醫療器械
ZH.	医疗器械
CS.	Zdravotnický prostředek
KK.	Медициналық құрылғы
PT.	Dispositivo médico

PT. Não contém látex
RO. Fără latex
AR. خالي من اللاتكس
FI. Ei sisällää luonnonkumilateksia

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DESCRIPTION

The “FLOW REVERSE” (code R9900101) manufactured by RAND is a disposable accessory for use in Hyperthermic Perfusion procedures and it is particularly useful to improve the fluid circulation and the uniform temperature distribution inside the peritoneal cavity.

Not made with natural rubber latex. Made of non-pyrogenic, non-toxic and biocompatible materials according to ISO 10993 standards and in relation with the intended use.

TECHNICAL FEATURES

Recommended perfusion flow: up to 2000 ml/min

Operating temperature: 37 to 46°C

Overall priming volume: about 40 ml

INTENDED USE

Disposable accessory intended to reverse the Patient's inlet and outlet flows during Hyperthermic Perfusion procedures.

SIDE EFFECTS AND CONTRAINDICATIONS

The device is contraindicated for procedures outside of the intended use.

There is no evidence of side effects specifically related to the use of the product.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer user's manual.

SAFETY INFORMATION

The information to draw the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device are indicated in the text as follows:

WARNING

Indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur both during proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such case.

PRECAUTION

Indicates any special precaution than the user must exercise for safe and effective use of the device.

For further information and/or in case of complaints, contact Rand directly or the authorised local representative.

WARNING

- The User must read and understand all User information prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- The device is intended for professional use only. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.
- Visually inspect and carefully check the devices before use. Transport and/or storage conditions not conform to the specifications may have caused damage to the devices.

- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the devices if sterility cannot be guaranteed.
- Check the expiry date shown on the label attached. Do not use the devices after this date.
- The devices must be used immediately after opening the sterile packaging.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo any further processing. Resterilization may not be effective and may compromise the integrity of the device thus resulting in unintended injury.
- Tubing should be connected in such a manner as to prevent kinks or narrowings.
- Avoid alcohol, anaesthetic fluids, or corrosive solvents from coming into contact with the devices as they may jeopardize the structural integrity.

PRECAUTION


- After use, dispose of the device in accordance with the applicable regulations in force in the country where it is used.

INSTALLATION

WARNING


Perform the procedure using an aseptic technique.

1. Connect the quick connectors of Flow Reverse to the Hang&Go or HIPP (F) equipment lines (precisely, to the CL2 external line and to the PM2 inlet line) on one side, and to the Inlet/Outlet table pack lines on the other side.

Note: the Flow Reverse does not have a preferential way of connection.


2. The starting configuration (normal flow) is with:
WHITE CLAMPS: OPEN
RED CLAMPS: CLOSED

3. To reverse the flow do the following steps in sequence:
 - i. Open the RED clamps
 - ii. Close the WHITE clamps

OPEN the RED clamps

■ ENGLISH ■

WARNING

- In order to avoid pressure alarms triggering on the Performer equipment, anytime the flow is reverted always **OPEN** the clamps **BEFORE** closing the other two.
- Close (or open) only clamps of the same color at the same time.
- Anytime the flow is reverted (i.e. anytime the User changes the configuration of the Flow Reverse from normal to reversed or viceversa), T1 and T2 terminals at the HTS temperature module of PERFORMER equipment must be exchanged too.

WARNING

The above mentioned devices have been designed for use up to a maximum pressure of 500 mmHg; Rand is not responsible for improper setting of the parameters or accidental over-pressurization of the devices exceeding this safety limit.

CIRCUIT PRIMING AND TREATMENT PHASES

For all instructions and warnings for the correct execution of these phases, follow the indications provided in the Hang&Go kit or HIPP (F) leaflet and in the PERFORMER user's manual.

COMPONENTS REPLACEMENT

As a precaution a spare device should always be available during priming and therapy phases. If particular conditions (e.g. leaks) potentially affecting the efficacy of the therapy or the patient/user's safety are detected on the Flow Reverse, interrupt the treatment and proceed to its replacement as follows:

- i. Activate the BY-PASS on the Performer equipment (refer to the user's manual of the equipment);
- ii. Clamp the four tubes connected to the Flow Reverse as close as possible to the quick connectors, using either the provided manual clamps or Klemmer clamps;
- iii. Close the four clamps of the Flow Reverse;
- iv. Disconnect the four quick connectors;
- v. Replace the device as described in point 1 of the section "Installation" above;
- vi. Open the clamps of the connected tubes;
- vii. Open the two clamps of the Flow Reverse that restore the configuration prior to the replacement;
- viii. Reactivate the "CIRCULATION" phase (refer to the Performer user's manual).

MEDICAL DEVICES TO BE USED IN COMBINATION

The FLOW REVERSE manufactured by RAND is designed to be used in connection with the tubing of the "HANG&GO HT" kit (code R9900067), "HANG&GO HT BASIC" (code R9900088 and R9900090) and "HIPP (F)" Line set (code R9900033).

WARNING

The compatibility of use of the devices listed in this instruction leaflet has been appropriately proved and verified through in vitro and in vivo tests. Rand will not be liable for any

consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to Rand or to its distributor and the competent authority of the Member State in which the user is established.

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RAND warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RAND warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RAND cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RAND, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RAND undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RAND up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RAND or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RAND disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RAND, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.