*Instructions for use
Istruzioni per l'uso
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Gebruiksaanwijzing
Használati utasítás
Bruksanvisning
Инструкция по использованию
Bruksanvisning
Kullanım talimatları*

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DESCRIPTION

The "FLOW REVERSER" (code R9900136) manufactured by RAND is a disposable accessory for use in Hyperthermic Perfusion procedures and it is particularly useful to improve the fluid circulation and the uniform temperature distribution inside the peritoneal cavity.

Not made with natural rubber latex. Made of non-pyrogenic, non-toxic and biocompatible materials according to ISO 10993 standards and in relation with the intended use.

INTENDED USE

Disposable accessory intended to invert the Patient's inlet and outlet flows during Hyperthermic Perfusion procedures.

SIDE EFFECTS AND CONTRAINDICATIONS

The device is contraindicated for procedures outside of the intended use. There is no evidence of side effects specifically related to the use of the product.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer 3 user's manual.

SAFETY INFORMATION

The information to draw the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device are indicated in the text as follows:

WARNING: indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur both during proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such case.

PRECAUTION: indicates any special precaution than the user must exercise for safe and effective use of the device.

For further information and/or in case of complaints, contact Rand directly or the authorised local representative.

WARNING

- The User must read and understand all User information prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- The device is intended for professional use only. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.
- Visually inspect and carefully check the devices before use. Transport and/or storage conditions not conform to the specifications may have caused damage to the devices.
- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the devices if sterility cannot be guaranteed.
- Check the expiry date shown on the label attached. Do not use the devices after this date.
- The devices must be used immediately after opening the sterile packaging.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo any further processing. Resterilization may not be effective and may compromise the integrity of the device thus resulting in unintended injury.
- Tubing should be connected in such a manner as to prevent kinks or narrowings.
- Avoid alcohol, anaesthetic fluids, or corrosive solvents from coming into contact with the devices as they may jeopardize the structural integrity.

PRECAUTION


- After use, dispose of the devices in accordance with the applicable regulations in force in the country where the device is used for biologically hazardous waste.

INSTALLATION


WARNING

Perform the procedure using an aseptic technique.

1. Connect the quick connectors of Flow Reverser to the Hang&Go 3 equipment lines, precisely to the CL3 lines and to the PM2 inlet line on one side, and to the Inlet/Outlet table pack lines on the other side.

2. The starting configuration (normal flow) is with:
WHITE CLAMPS: OPEN
RED CLAMPS: CLOSED

INVERTING THE FLOW

Refer to the chapter "FLOW REVERSER" in the Performer 3 user's manual.

WARNING

- In order to avoid pressure alarms triggering on the Performer equipment, anytime the flow is inverted always OPEN the clamps BEFORE closing the other two.
- Close (or open) only clamps of the same color at the same time.

CIRCUIT PRIMING AND TREATMENT PHASES

For all instructions and warnings for the correct execution of the treatment phases, follow the indications provided in the Hang&Go 3 leaflet and in the Performer 3 user's manual.

PRIMING VOLUME

40 mL.

COMPONENTS REPLACEMENT

As a precaution a spare device should always be available during priming and therapy phases. If particular conditions (e.g. leaks) potentially affecting the efficacy of the therapy or the patient/user's safety are detected on the Flow Reverser, interrupt the treatment and proceed to its replacement as follows:

- i. Activate the BY-PASS on the Performer equipment (refer to the user's manual of the equipment);
- ii. Clamp the 6 tubes connected to the Flow Reverser as close as possible to the quick connectors, using either the provided manual clamps or Klemmer clamps;
- iii. Close the 4 clamps of the Flow Reverser;
- iv. Disconnect the 6 quick connectors;
- v. Replace the device as described in point 1 of the section "Installation" above;
- vi. Open the clamps of the connected tubes;
- vii. Open the 2 clamps of the Flow Reverser that restore the configuration prior to the replacement;
- viii. Reactivate the "CIRCULATION" phase (refer to the Performer user's manual).

MEDICAL DEVICES TO BE USED IN COMBINATION

The Flow Reverser is designed to be used in connection with the tubing of the Hang&Go 3 (code R9900141).

WARNING

Rand will not be liable for any consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to Rand or to its distributor and the competent authority of the Member State in which the user is established.

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RAND warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RAND warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RAND cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RAND, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RAND undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RAND up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RAND or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RAND disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RAND, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.