Instructions for use
Istruzioni per l'uso
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Gebruiksaanwijzing
Használati utasítás
Brugsanvisninger
Инструкция по использованию
Bruksanvisning
Kullanım talimatları
使用说明书

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DESCRIPTION

Sterile single-use devices designed and intended to be used in combination with RanD Performer 3 equipment to enable the entire system to reach the intended purpose, i.e. to provide warmed fluids to the patient in the Hyperthermic Intra-Peritoneal / Intra-pleural Chemotherapy (also called HIPEC) in the oncology field. Specifically, they are extracorporeal circulation tubing sets realized with biocompatible medical grade plastic materials intended to physically connect the patient to the dedicated equipment in order to allow the circulation of the warmed fluids from the equipment to the patient and back, in a continuous closed loop circulation.

Hang&Go 3 (code R9900141) is composed of:

- 1 Hang&Go 3 Equipment (code R9900150), line set for hyperthermic perfusion
- 1 Hang&Go 3 Table Pack (code R9900151), line set for hyperthermic perfusion

Hang&Go 3 Plus (code R9900140) is a procedure pack produced by RAND and is composed of:

- 1 Hang&Go 3 Equipment (code R9900150), line set for hyperthermic perfusion
- 1 Hang&Go 3 Table Pack (code R9900151), line set for hyperthermic perfusion
- 1 Flow Reverser (code R9900136)
- 3 Ch24 round silicone catheters (code R9900071)
- 3 Redax round coaxial drains (code 24727)
- 3 Metko disposable temperature probes (code FMT400/AOR-D5)

Hang&Go 3 DR (code R9900160) is a procedure pack produced by RAND and is composed of:

- 1 Hang&Go 3 Equipment (code R9900150), line set for hyperthermic perfusion
- 1 Hang&Go 3 Table Pack (code R9900151), line set for hyperthermic perfusion
- 3 Ch24 round silicone catheters (code R9900071)
- 3 Redax round coaxial drains (code 24727)

Hang&Go 3 HV (code R9900143) is composed of:

- 1 Hang&Go 3 Equipment (code R9900150), line set for hyperthermic perfusion
- 1 Hang&Go 3 Table Pack (code R9900151), line set for hyperthermic perfusion
- 1 10000mL Waste Bags (code R9900111)
- 6 Table pack 1/4 connectors (code R9900122)

INTENDED USE

The Hang&Go 3, the Hang&Go 3 Plus, the Hang&Go 3 DR and the Hang&Go 3 HV are intended to recirculate, filtrate and perfuse clear solutions in the thoracic or abdominal cavity. They are intended for use up to 6 hours, in regional cancer hyperthermic therapies following surgical reduction.

SIDE EFFECTS AND CONTRAINDICATIONS

There is no evidence of side effects specifically related to the use of the product. The device is contraindicated for procedures outside of the intended use.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer 3 user's manual.

SAFETY INFORMATION

The information to draw the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device are indicated in the text as follows:

WARNING: indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur both during proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such case.

PRECAUTION: indicates any special precaution than the user must exercise for safe and effective use of the device

For further information and/or in case of complaints, contact RanD directly or the authorised local representative.

WARNING

- The user must read and understand all user information prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- The devices are intended for professional use only. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.
- Visually inspect and carefully check the devices before use. Transport and/or storage conditions non conform to the specifications may have caused damage to the devices.
- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the devices if sterility cannot be guaranteed.
- Check the expiry date shown on the label attached. Do not use the devices after this date.
- The devices must be used immediately after opening the sterile packaging.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo any further processing. Resterilization may not be effective and may compromise the integrity of the device thus resulting in unintended injury.
- Tubing should be connected in such a manner as to prevent kinks or narrowings.
- Avoid alcohol, anaesthetic fluids, or corrosive solvents from coming into contact with the devices as they may jeopardize the structural integrity

PRECAUTION

- As a precaution, a back-up Hang&Go 3 should always be available during priming and therapy phases. If you believe that a particular event (such as insufficient performances, leaks, filter saturation, etc.) may affect the patient/user safety, interrupt the treatment as explained in the Performer 3 user's manual (Section: "How to force interruption of the treatment during an HIPEC procedure") and proceed to the device replacement.
- Never unpack or prepare the table pack lines close to the patient's surgical wound. Detachable parts not well secured could fall into the patient's cavity without being noticed.
- After use, dispose of the devices in accordance with the applicable regulations in force in the country where the device is used for biologically hazardous waste.

INSTALLATION

Refer to "Tubing set-up" in the Performer 3 user's manual.

WARNING

Perform the procedure using an aseptic technique.

CIRCUIT PRIMING, TREATMENT EXECUTION AND ENDING

Refer to the following chapters in the Performer 3 user's manual:

- PREPARATION PHASE
- PATIENT FILLING PHASE
- CIRCULATION PHASE
- EMPTYING PHASE
- RINSING PHASE
- END OF THE PROCEDURE

PRIMING VOLUME

1.5 L of sterile solution (including 1 L to prime reservoir and lines and minimum 0.5 L for Patient volume).

PRECAUTION

During perfusion it is advisable to check the state of saturation of the filters integrated in the semi-rigid reservoir hung on the equipment scale. Irrespectively on the patient's pathology, treatment duration and surgical approaches (closed or open abdomen perfusion), make sure that the filters are pervious throughout the treatment phase by visually cross-checking.

PRECAUTION

RanD advises to reach the right treatment parameters gradually, to ensure a proper therapy start as well as a correct treatment execution.

TREATMENT, DRUG INJECTION (OPTIONAL):

RanD suggests injecting medication through the (blue) Clave connectors integrated at different points in the circuit.

In order to obtain the best therapeutic result, it is advisable to inject drug only when the perfused district has already reached the required treatment temperature.

WARNING

The compatibility of use of the device with some commonly used medications - at the discretion of the physician - in intraperitoneal hyperthermic perfusion treatments has been appropriately verified by careful review of scientific literature and analysis of clinical test data relating to use of commercial medications based on the following active ingredients of consolidated use: CARBOPLATIN, CISPLATIN, MITOMYCIN C, ADRIAMYCIN/DOXORUBICIN, IRINOTECAN, MELPHALAN and OXALIPLATIN. RanD does not assume any responsibility for use of the devices in combination with medications or preparations not previously evaluated.

MEDICAL DEVICES TO BE USED IN COMBINATION

The Hang&Go 3, the Hang&Go 3 Plus and the Hang&Go 3 HV are designed for use in combination with the RanD Performer 3 equipment.

Ch 24 round silicone Catheters (code R9900071) have been specifically designed for use in combination with the "Patient Inlet" (code R2500286) and "Patient Outlet" (code R2500285) lines of the Hang&Go 3 Table Pack.

The Flow Reverser (code R9900136) is an accessory intended to be connected between the Hang&Go 3 Equipment lines and the Hang&Go 3 Table Pack lines. For set-up and use of this accessory, please refer to the instructions for use provided in the Flow Reverser package.

The Redax round coaxial drains (code 24727) are intended to be connected to the withdrawal lines of the Hang&Go 3 table pack. They are not intended for use in infusion.

In case the volume to recover from patient's body cavity is greater than 8 liters (maximum capacity of the reservoir bag), connect first the 5 liters waste bag and then the 10 liters waste bag (included in the Hang&Go 3 HV kit).

WARNING

RanD will not be liable for any consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to RanD or to its distributor and the competent authority of the Member State in which the user is established

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RanD warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RanD warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with

■ ENGLISH ■

them by a qualified user and before any expiry date indicated on the packaging. However, RanD cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RanD, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RanD undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RanD up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RanD or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RanD disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RanD, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.