*Instructions for use
Istruzioni per l'uso
Mode d'emploi
Gebrauchsanweisung
Gebruiksaanwijzing
Brugsanvisninger*

R9900119 | HANG&GO BASIC HV

EN	<i>Kit for Hyperthermic Perfusion</i>
IT	<i>Kit per Perfusione Ipertermica</i>
FR	<i>Kit Pour Perfusion Hyperthermique</i>
DE	<i>Set für hypertherme Perfusion</i>
NL	<i>Kit voor hyperthermische perfusie</i>
DA	<i>Sæt til hypertermisk perfusion</i>

EN. Do not reuse
IT. Non riutilizzare
FR. Ne pas réutiliser
DE. Nicht zur Wiederverwendung
ES. No reutilizar
TR. Tek kullanımlık
NL. Niet opnieuw gebruiken
RU. Обратитесь к инструкции по применению
PL. Do jednorazowego użytku
SR. Nije za višekratnu upotrebu
LT. Nenaudokite pakartotinai
SV. Får inte återanvändas
BG. Да не се използва повторно
DA. Må ikke genbruges
HR. Nije za višekratnu upotrebu
ZH. 请不要重复使用
ZH. 不可重複使用
CS. Nepoužívejte opakovaně
KK. Қайта пайдалануға болмайды
PT. Não reutilizar
RO. A nu se reutiliza
AR. لا تُعد استخداماً

EN. Content
IT. Contenuto
FR. Contenu
DE. Inhalt
ES. Contenido
TR. İçerik
NL. Inhoud
RU. Содержит
PL. Zawiera
SR. Sadržaji
LT. Sudėtis
SV. Innehåll
BG. Съдържание
DA. Indhold
HR. Sadržaj
ZH. 内容
ZH. 潤德
CS. Obsah
KK. Содержимое
PT. Conteúdo
RO. Conținut
AR. المحتوى

EN. Non-pyrogenic
IT. Apirogeno
FR. Apyrogène
DE. Pyrogenfrei
ES. Apirogeno
TR. Apirojen
NL. Niet-pyrogeen
RU. Апиrogenно
PL. Niepirogeny
SR. Ne stvara toplotu
LT. Nepirogeniškas
SV. Pyrogenfri
BG. Апиrogenно
DA. Ikke-pyrogen
HR. Nepirogeno
ZH. 无热源的
ZH. 無致熱原
CS. Nepyrogní
KK. Апиrogenді
PT. Não pirogénico
RO. Apirogen
AR. غير بيريوجيني

EN. Quantity
IT. Quantità
FR. Quantité
DE. Menge
ES. Cantidad
TR. Adet
NL. Aantal
RU. Количество
PL. Ilość
SR. Količina
LT. Kiekis
SV. Antal
BG. Количество
DA. Antal
HR. Količina
ZH. 数量
ZH. 數量
CS. Množství
KK. Саны
PT. Quantidade
RO. Cantitate
AR. الكمية

EN. Keep dry
IT. Teme l'umidità
FR. Craint l'humidité
DE. Trocken lagern
ES. Matener seco
TR. Kuru saklayın
NL. Droog bewaren
RU. Не допускать воздействия влаги
PL. Przechowywać w suchym miejscu
SR. Čuvati na suvom
LT. Laikykite sausai
SV. Förvaras torrt
BG. Да се държи на сухо място
DA. Skal holdes tørt
HR. Čuvati na suhom
ZH. 请保持干燥
ZH. 保持乾燥
CS. Uchovávejte v suchu
KK. Құрғақ сақтаңыз
PT. Manter seco
RO. A se păstra ferit de umezeală
AR. حافظ على جفافه

EN. Date of manufacture
IT. Data di fabbricazione
FR. Date de fabrication
DE. Herstellungsdatum
ES. Fecha de fabricación
TR. Üretim tarihi
NL. Fabricagedatum
RU. Дата производства
PL. Data produkcji
SR. Datum proizvodnje
LT. Pagaminimo data
SV. Tillverkningsdatum
BG. Дата на производство
DA. Fremstillingsdato
HR. Datum proizvodnje
ZH. 生产日期
ZH. 製造日期
CS. Datum výroby
KK. Өндіру күні
PT. Data de fabrico
RO. Data fabricației
AR. تاريخ التصنيع

EN. Do not re-sterilize
IT. Non ristilizzare
FR. Ne pas restériliser
DE. Nicht resterilisieren
ES. No volver a esterilizar
TR. Tekrar steril etmeyin
NL. Niet opnieuw steriliseren
RU. Не стерилизовать повторно
PL. Nie sterylizować
SR. Ne sterilizovati više puta
LT. Nesterilizuoti pakartotinai
SV. Får inte återsteriliseras
BG. Да не се стерилизира повторно
DA. Må ikke gensteriliseres
HR. Ponovna sterilizacija nije dopuštena
ZH. 请勿重复消毒
ZH. 不可重複滅菌
CS. Nesterilizujte opakovaně
KK. Қайта стерилдеуге болмайды
PT. Não reesterilizar
RO. A nu se resteriliza
AR. لا تُعد تعقيم

EN. Use by
IT. Utilizzare entro
FR. Utiliser jusqu'à
DE. Verwendbar bis
ES. Fecha de caducidad
TR. Son kullanım tarihi
NL. Te gebruiken vóór
RU. Использовать до
PL. Zużyć do
SR. Rok upotrebe
LT. Tinka naudoti iki
SV. Använd före
BG. Годен до
DA. Anvendes inden
HR. Upotrijebiti do
ZH. 使用者
ZH. 保存期限
CS. Spotřebujte do
KK. Жарамдылық мерзімі
PT. Data de validade
RO. A se utiliza înainte de
AR. يستخدم من قبل

EN. Batch code
IT. Numero di lotto
FR. Numéro de lot
DE. Chargenbezeichnung
ES. Número de lote
TR. Lot numarası
NL. Partijnummer
RU. Номер партии
PL. Kod partii
SR. Šifra serije
LT. Partijos numeris
SV. Batchnummer
BG. Партиден номер
DA. Partikode
HR. Šifra serije
ZH. 生产批号
ZH. 製造批號
CS. Kód šarže
KK. Жинақ коды
PT. Código de lote
RO. Număr lot
AR. رمز التشغيل

EN. Manufacturer
IT. Fabricante
FR. Fabricant
DE. Hersteller
ES. Fabricante
TR. Üretici
NL. Fabrikant
RU. Производитель
PL. Producent
SR. Proizvođač
LT. Gamintojas
SV. Tillverkare
BG. Производител
DA. Producent
HR. Proizvođač
ZH. 生产厂家
ZH. 製造廠
CS. Výrobce
KK. Өндіруші
PT. Fabricante
RO. Producător
AR. الشركة المصنعة

EN. Assembled by
IT. Assemblato da
FR. Assemblé par
DE. Zusammengefügt von
ES. Ensamblado por
TR. Birleştirici
NL. Geassembleerd door
RU. Упаковщик
PL. Złożono
SR. Sastavio
LT. Surinko
SV. Monterad av
BG. Сглобено от
DA. Samlet af
HR. Sastavio
ZH. 组装
CS. Sestavil
KK. Құрастырушы
PT. Montado por
RO. Asamblat de
AR. التجميع من قبل

EN. Catalogue number
IT. Numero di catalogo
FR. Référence du catalogue
DE. Bestellnummer
ES. Número de catálogo
TR. Referans numarası
NL. Catalogusnummer
RU. Номер по каталогу
PL. Numer katalogowy
SR. Kataloški broj
LT. Katalogo numeris
SV. Katalognummer
BG. Каталоген номер
DA. Katalognummer
HR. Broj u katalogu
ZH. 目录号
ZH. 產品型號
CS. Katalogové číslo
KK. Каталог нөмірі
PT. Número de catálogo
RO. Număr de catalog
AR. رقم الكatalog

EN. Keep away from heat
IT. Tenere lontano da sorgenti di calore
FR. Craint la chaleur
DE. Vor Wärmequellen fernhalten
ES. Mantener lejos de fuentes de calor
TR. Isı kaynaklarından uzak tutun
NL. Uit de buurt van warmtebronnen houden
RU. Не допускать воздействия тепла
PL. Przechowywać z dala od źródeł ciepła
SR. Držati dalje od toplote
LT. Saugokite nuo šilumos
SV. Skydda från värme
BG. Да се пази от топлина
DA. Må ikke udsættes for varme
HR. Držati dalje od topline
ZH. 远离高温
ZH. 避免日照
CS. Chraňte před teplem
KK. Ыстықтан алыс ұстаңыз
PT. Manter afastado do calor
RO. A se păstra departe de sursele de căldură
AR. أبقيه بعيداً عن الحرارة

STERILE EO

EN. Sterile - Ethylene oxide sterilized
IT. Sterile - Sterilizzato ad ossido d'etilene
FR. Stérile - Stérilisé à l'oxyde d'éthylène
DE. Steril - Mit Äthylenoxid sterilisiert
ES. Estéril - Esterilizado por óxido de etileno
TR. Steril - Etilen Oksit
NL. Steriel - Gesteriliseerd met ethylenoxide
RU. Стерильно - Стерилизовано окисью этилена
PL. Sterylne - Wysterylizowane tlenkiem etylenu
SR. Sterilno - Sterilisano etilen oksidom
LT. Sterilu - Sterilizuota etileno oksidu
SV. Steril - Steriliserad med etylenoxid
BG. Стерильно - Стерилизирано с этиленов оксид
DA. Steril - Steriliseret med ethylenoxid
HR. Sterilno - Sterilizirano etilen-oksidiom
ZH. 消毒 - 环氧乙烷消毒法
ZH. EO滅菌
CS. Sterilní - Sterilizováno ethylenoxidem
KK. Стерильді - Этилен тотығымен стерильденген
PT. Estéril - Esterilizado por óxido de etileno
RO. Steril - Sterilizat cu oxid de etilenă
AR. معقم - معقم بأكسيد الإيثيلين

EN. Not made with natural rubber latex
IT. Non fatto con lattice di gomma naturale
FR. Latex exempt
DE. Latex-frei
ES. Libre de látex
TR. Lateks içermez
NL. Latexvrij
RU. Не содержит натуральный каучуковый латекс
PL. Nie zawiera lateksu
SR. Ne sadrži lateks
LT. Be latekso
SV. Latexfri
BG. Не съдържа латекс
DA. Latex-fri
HR. Bez lateksa
ZH. 请与乳胶隔离
ZH. 非乳膠製品
CS. Neobsahuje latex
KK. Латексиз
PT. Não contém látex
RO. Fără latex
AR. خالي من اللاتكس

EN. Temperature limitation
IT. Limiti di temperatura
FR. Limites de température
DE. Zulässiger Temperaturbereich
ES. Límite de temperatura
TR. Sıcaklık limitleri
NL. Temperatuurlimieten
RU. Диапазон температур
PL. Ograniczenia temperatury
SR. Temperaturno ograničenje
LT. Temperatūros apribojimas
SV. Temperaturbegränsning
BG. Температурно ограничение
DA. Temperaturbegrænsning
HR. Ograničenje temperature
ZH. 温度限制
ZH. 溫度限制
CS. Teplotní omezení
KK. Температур шери
PT. Limite de temperatura
RO. Limita de temperatură
AR. حدود درجة الحرارة

EN. Caution, consult accompanying documents
IT. Attenzione, vedere le istruzioni per l'uso
FR. Attention voir notice d'instructions
DE. Achtung, Begleiddokumente beachten
ES. Atención, ver instrucciones de uso
TR. Dikkat, kullanım talimatlarına bakın
NL. Attentie, de bijgevoegde documenten raadplegen
RU. Обратитесь к инструкции по применению
PL. Uwaga! Zapoznać się z dołączonymi dokumentami
SR. Pažnja, pogledajte prateću dokumentaciju
LT. Dėmesio, perskaitykite pridedamus dokumentus
SV. Varning – läs de medföljande dokumenten
BG. Внимание, направете справка в придружаващата документация
DA. Forsigtig, rådfør dig med de medfølgende dokumenter
HR. Oprez, pogledati prateće dokumente
ZH. 注意. 请阅读所附文件
ZH. 注意事項，請參考隨附文件
CS. Upozornění: prostudujte si průvodní dokumentaci
KK. Ескерту, бірге келген құжаттарды қараңыз
PT. Cuidado, consultar os documentos fornecidos
RO. Atenție, consultați documentele însoțitoare
AR. احترس، راجع الوثائق المرفقة

EN. Do not use if package is damaged
IT. Non utilizzare se la confezione è danneggiata
FR. Ne pas utiliser si l'emballage est endommagé
DE. Inhalt beschädigter Packung nicht verwenden
ES. No usar si el paquete está dañado
TR. Paket hasar görmüşse kullanmayın
NL. Niet gebruiken als de verpakking beschadigd is
RU. Не использовать при повреждении упаковки
PL. Nie używać, jeśli opakowanie jest uszkodzone
SR. Ne koristiti ukoliko je pakovanje oštećeno
LT. Nenaudokite, jeigu pažeista pakuotė
SV. Använd inte om förpackningen är skadad
BG. Да не се използва, ако опаковката е повредена
DA. Må ikke anvendes, hvis emballagen er beskadiget
HR. Ne upotrebljavati ako je pakiranje oštećeno
ZH. 包装出现损坏的情况下请不要使用
ZH. 包裝破損不可使用
CS. Nepoužívejte, pokud je poškozený obal
KK. Орам зақымдалған болса, пайдаланбаңыз
PT. Não utilizar se a embalagem estiver danificada
RO. A nu se utiliza dacă ambalajul este deteriorat
AR. لا تستخدمه في حالة تلف الحزمة

EN. Medical device
IT. Dispositivo medico
FR. Dispositif médical
DE. Medizinprodukt
ES. Producto sanitario
TR. Tıbbi cihaz
NL. Medisch hulpmiddel
RU. Медицинское изделие
PL. Wyrób medyczny
SR. Medicinsko sredstvo
LT. Medicinos priemonė
SV. Medicinteknisk produkt
BG. Медицинско изделие
DA. Medicinsk udstyr
HR. Medicinski proizvod
ZH. 醫療器械
ZH. 医疗器械
CS. Zdravotnický prostředek
KK. Медициналық құрылғы
PT. Dispositivo médico
RO. Dispozitiv medical
AR. جهاز طبي

EN. Single sterile barrier system
IT. Sistema a singola barriera sterile
FR. Système à simple barrière stérile
DE. Einfaches Sterilbarrieresystem
ES. Sistema de única barrera estéril
TR. Tek steril bariyer sistemi
NL. Enkel steriel barrièresysteem
RU. Одинарная стерильная барьерная система
PL. System pojedynczej bariery sterylnej
SR. Sistem sa jednom sterilnom barijerom
LT. Vienguba sterilių barjerų sistema
SV. Enkelt sterilt barriärsystem
BG. Единична стерилна бариерна система
DA. Enkelt sterilt barrieresystem
HR. Sustav jednostruke sterilne obloge
ZH. 單層無菌屏障系統
ZH. 单层无菌屏障系统
CS. Systém s jednou sterilní bariérou
KK. Бир стерильді тосқауыл жүйесі
PT. Sistema de barreira estéril único
RO. Sistem de barieră sterilă simplă
AR. نظام عازل معقم ومنفرد

EN. Double sterile barrier system
IT. Sistema a doppia barriera sterile
FR. Système à double barrière stérile
DE. Doppeltes Sterilbarrieresystem
ES. Sistema de doble barrera estéril
TR. Çift steril bariyer sistemi
NL. Dubbel steriel barrièresysteem
RU. Двойная стерильная барьерная система
PL. System podwójnej bariery sterylnej
SR. Sistem sa dve sterilne barijere
LT. Dviguba sterilių barjerų sistema
SV. Dubbelt sterilt barriärsystem
BG. Двойна стерилна бариерна система
DA. Dobbelt sterilt barrieresystem
HR. Sustav dvostruke sterilne obloge
ZH. 雙層無菌屏障系統
ZH. 双层无菌屏障系统
CS. Systém se dvěma sterilními bariérami
KK. Екі стерильді тосқауыл жүйесі
PT. Sistema de barreira estéril duplo
RO. Sistem de barieră sterilă dublă
AR. نظام عازل معقم مزدوج

EN. Single sterile barrier system with protective packaging inside
IT. Sistema a singola barriera sterile con confezionamento protettivo interno
FR. Système à simple barrière stérile avec emballage de protection à l'intérieur
DE. Einfaches Sterilbarrieresystem mit Schutzverpackung innen
ES. Sistema de única barrera estéril con embalaje protector en el interior
TR. İçinde koruyucu ambalaj bulunan tek steril bariyer sistemi
NL. Enkel steriel barrièresysteem met beschermende verpakking aan de binnenzijde
RU. Одинарная стерильная барьерная система с защитной упаковкой внутри
PL. System pojedynczej bariery sterylnej z ochronnym opakowaniem wewnętrznym
SR. Sistem sa jednom sterilnom barijerom i unutrašnjim zaštitnim pakovanjem
LT. Vienguba sterilaus barjero sistema su apsaugine pakuote viduje
SV. Enkelt sterilt barriärsystem med skyddande förpackningar inuti
BG. Единична стерилна бариерна система с вътрешна защитна опаковка
DA. Enkelt sterilt barrieresystem med indvendig beskyttende emballage
HR. Sustav jednostruke sterilne obloge s unutarnjom zaštitnom ambalažom
ZH. 單層無菌屏障系統，內有保護性包裝
ZH. 单层无菌屏障系统及内层防护包装
CS. Systém s jednou sterilní bariérou s ochranným obalem uvnitř
KK. Ішінде қорғаныс орамы бар бір стерильді тосқауыл жүйесі
PT. Sistema de barreira estéril único com embalagem protetora interior
RO. Sistem de barieră sterilă simplă cu ambalaj de protecție în interior
AR. نظام عازل معقم ومنفرد بعنوة واقية من الداخل

EN. Single sterile barrier system with protective packaging outside
IT. Sistema a singola barriera sterile con confezionamento protettivo esterno
FR. Système à simple barrière stérile avec emballage de protection à l'extérieur
DE. Einfaches Sterilbarrieresystem mit Schutzverpackung außen
ES. Sistema de única barrera estéril con embalaje protector en el exterior
TR. Dışında koruyucu ambalajlı tek steril bariyer sistemi
NL. Enkel steriel barrièresysteem met beschermende verpakking aan de buitenzijde
RU. Одинарная стерильная барьерная система с защитной упаковкой снаружи
PL. System pojedynczej bariery sterylnej z ochronnym opakowaniem zewnętrznym
SR. Sistem sa jednom sterilnom barijerom i spoljašnjim zaštitnim pakovanjem
LT. Vienguba sterilaus barjero sistema su apsaugine pakuote išorėje
SV. Enkelt sterilt barriärsystem med skyddande förpackning utanför
BG. Единична стерилна бариерна система с външна защитна опаковка
DA. Enkelt sterilt barrieresystem med udvendig beskyttende emballage
HR. Sustav jednostruke sterilne obloge s vanjskom zaštitnom ambalažom
ZH. 單層無菌屏障系統，外部有保護性包裝
ZH. 单层无菌屏障系统及外层防护包装
CS. Systém s jednou sterilní bariérou s ochranným obalem vně
KK. Сыртында қорғаныс орамы бар бір стерильді тосқауыл жүйесі
PT. Sistema de barreira estéril único com embalagem protetora exterior
RO. Sistem de barieră sterilă simplă cu ambalaj de protecție în exterior
AR. نظام عازل معقم ومنفرد بعنوة واقية من الخارج

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■ ENGLISH ■

DESCRIPTION

Sterile single-use devices designed and intended to be used in combination with Rand Performer HT equipment to enable the entire system to reach the intended purpose, i.e. to provide warmed fluids to the patient in the Hyperthermic Intra-Peritoneal / Intra-pleural Chemotherapy (also called HIPEC) in the oncology field.

Specifically, they are extracorporeal circulation tubing sets realized with biocompatible medical grade plastic materials intended to physically connect the patient to the dedicated equipment in order to allow the circulation of the warmed fluids from the equipment to the patient and back, in a continuous closed loop circulation.

Each "HANG&GO BASIC HV" (code R9900119) contains the following devices:

- 1 "Hang&Go HT Equipment" line set (code R9900146),
- 1 "Hang&Go Table Pack" line set (code R9900149),
- 1 "10000mL Waste Bags" (code R9900111).

TECHNICAL FEATURES

Recommended perfusion flow: up to 2000 ml/min

Operating temperature: 37 to 46°C

Overall priming volume of the extracorporeal circuit (including the table pack lines): 500 ml

Overall priming volume: 1000 ml (including a minimum 500 ml level within the filtered reservoir).

The HANG&GO BASIC HV must be used on Performer HT with software version 3.0.5.E or greater, featuring the "high volume" option.

INTENDED USE

Disposable devices intended to recirculate, filtrate and perfuse clear solutions, also enriched with chemotherapeutic agents, in the thoracic or abdominal cavity. The devices are intended for use up to 6 hours, in regional cancer hyperthermic therapies following surgical reduction.

SIDE EFFECTS AND CONTRAINDICATIONS

There is no evidence of side effects specifically related to the use of the product. The device is contraindicated for procedures outside of the intended use.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer equipment user's manual.

SAFETY INFORMATION

The information to draw the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device are indicated in the text as follows:

WARNING

Indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur both during proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such case.

PRECAUTION

Indicates any special precaution than the user must exercise for safe and effective use of the device.

Below follows some general safety information for proper use of the device. Furthermore, specific safety-related information is given in the instructions for use where they may condition the operation to be performed.

For further information and/or in case of complaints, contact Rand directly or the authorised local representative.

WARNING

- The user must read and understand all user information prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- The devices are intended for professional use only. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.
- Visually inspect and carefully check the devices before use. Transport and/or storage conditions not conform to the specifications may have caused damage to the devices.
- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the devices if sterility cannot be guaranteed.
- Check the expiry date shown on the label attached. Do not use the devices after this date.
- The devices must be used immediately after opening the sterile packaging.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo any further processing. Resterilization may not be effective and may compromise the integrity of the device thus resulting in unintended injury.

■ ENGLISH ■

- Tubing should be connected in such a manner as to prevent kinks or narrowings.
- Avoid alcohol, anaesthetic fluids, or corrosive solvents from coming into contact with the devices as they may jeopardize the structural integrity.

PRECAUTION



- Never unpack or prepare the table pack lines close to the patient's surgical wound. Detachable parts not well secured could fall into the patient's cavity without being noticed.
- After use, dispose of the devices in accordance with the applicable regulations in force in the country where the device is used.

INSTALLATION

WARNING

Perform the procedure using an aseptic technique.

1. Remove the table pack line set R9900149 to access the "Hang&Go HT Equipment" line set (code R9900146).
2. Open the "Hang&Go HT Equipment" line set package and extract the "Hang&Go HT Reservoir" (code R2500273) from the sterile wrapping; check the position and integrity of each component.
3. Verify that the equipment is in the Home Screen.
4. Hook the filtrate reservoir onto the scale by means of the two rings on the head of the rigid part.
5. Close the manual clamp at the bottom of the soft bag.
6. Remove the blue marked line from its slot on the reservoir right side.
7. Carefully place the pump segment tubing into PM1 and PM2 pumps, paying attention to the fitting direction: bear in mind that all the pumps in the PERFORMER HT equipment turn clockwise. In order to facilitate the correct set up of the lines, each pump segment is marked with a yellow adhesive collar on the inlet side, which indicates the inlet side and the corresponding pump.
8. Insert the blue-marked soft tubing segment into the outer channel of the two-way clamp (CL1). The blue label must be positioned below the clamp.
9. Insert the tubing in the air sensor (SA1), paying attention to the correct closure of the sensor cover.
10. Insert the red-marked soft tubing segment into the outer channel of the two-way clamp (CL2). The red label must be positioned below the clamp.
11. Check the position and integrity of the membrane transducers used on the pressure monitoring lines; depending on use (negative or positive pressure monitoring), the deformable membrane may be on either side of the transducer. To facilitate the setting up procedure, the body of each transducer is clearly marked with a yellow label placed on the side where the membrane must be positioned for correct use. In addition, the label shows the name of the pressure gauge to which the pressure transducer must be connected (PR1, PR2, PR3 and PR5).
The following pictures show the correct positioning of the deformable membrane for monitoring **negative** pressure (PR1 and PR3) or **positive** pressure (PR2 and PR5):

If necessary, connect a syringe to the fluid isolator membrane connector and then inject or remove air until the membrane is correctly oriented.

12. Connect the pressure lines to the corresponding transducers on the machine.
13. Connect the male quick-connector of CL2 external line with the female quick-terminal at PM2 inlet.
14. Remove the two lines from their slots on the front of the reservoir rigid part.
15. Take out the heater bag (code R2500274) from "Hang&Go HT Equipment" line set, remove the sterile wrapping and place the heater bag in the heater compartment by means of the heater bag eyelets present at the four corners.
16. Carefully close the heater cover.
17. Connect the green male quick-connector at the heater bag inlet with the green female connector at the PM1 outlet.
18. Connect the brown male quick-connector at the heater bag outlet with the brown female quick-terminal of the line positioned on CL2 pinch-valve.
19. Enter Perfusion mode: the two clamps change their status.
20. Insert the two unmarked soft tubing segments into the applicable inner channels of CL1 and CL2 clamps.
21. Hang the perfusion solution bags on the IV pole.
22. Close the two manual clamps on the Priming line.
23. Insert the two spikes of the Priming line into the solution bags (if necessary, two additional spikes are provided in a bag stuck outside the "HANG&GO HT Reservoir" pouch, code R2500273).
24. Open the two manual clamps on the Priming line.
25. Verify proper setting-up of the circuit.

CIRCUIT PRIMING PROCEDURE

Carefully follow the step-by-step instructions provided in the PERFORMER HT user's manual for a correct priming of the disposable circuit.

PRECAUTION

- 1) For priming and rinsing procedures it is recommended to use at least 2000 ml of sterile saline solution.
- 2) During these phases, it is advisable to lightly tap the outside of the tubing in order to facilitate bleeding all the air from the circuit, especially from low flow segment.

THERAPY

PREPARATION, PATIENT FILLING PHASE:

After the applicable surgical preparation procedure is completed, including placement of infusion catheters, drainages and temperature probes in the patient's body cavity, the personnel operating in the sterile field shall connect:

- the "Infusion table pack" line code R2500237 to the infusion catheters
 - the "Withdrawal table pack" line code R2500236 to the withdrawal catheters
1. Disconnect the male quick-connector of CL2 external line from the female quick-connector at PM2 inlet.
 2. Connect the "Infusion table pack" line to the male quick-connector of CL2 external line.
 3. Connect the "Withdrawal table pack" line to the female quick-connector at PM2 inlet.
 4. After connecting the above mentioned two lines, unroll the integrated insulation sheaths over the withdrawal and infusion lines.
 5. Insert the temperature probe cable connectors into the corresponding HTS jacks (T1 through T8):
 - Connect the (red marked) temperature probe assembled in the "Infusion table pack" line to the T1 channel of the HTS module, by means of the appropriate cable;
 - Connect the (blue marked) temperature probe assembled in the "Withdrawal table pack" line to the T2 channel of the HTS module, by means of the appropriate cable;
 - Connect the remaining temperature probes to the corresponding channels of the HTS module (where appropriate);
 6. Make sure that the manual clamps on CL2 external line and on PM2 inlet line are open.

PRECAUTION

Select the appropriate volume and temperature on the equipment and load the volume of solution to perfuse the patient; in this phase, carefully ensure that the air coming from the withdrawal table pack is adequately held back and/or eliminated by the filtration system integrated in the reservoir hung on the device scale.

WARNING

The surgical team is directly responsible for the correct catheters placement.

WARNING

The above mentioned devices have been designed for use up to a maximum pressure of 500 mmHg; Rand is not responsible for improper setting of the parameters or accidental over-pressurization of the devices exceeding this safety limit.

TREATMENT, GENERAL NOTES:

- For a correct execution of this phase, follow the conditions and the indications provided with the PERFORMER HT user's manual.
- During perfusion it is advisable to check the state of saturation of the filters integrated in the semi-rigid reservoir hung on the equipment scale. Irrespectively on the patient's pathology, treatment duration and surgical approaches (closed or open abdomen perfusion), make sure that the filters are pervious throughout the treatment phase, by monitoring the visualized pressures on the PERFORMER HT screen and by visually cross-checking.

PRECAUTION

RAND advises to reach the right treatment parameters gradually, to ensure a proper therapy start as well as a correct treatment execution.

TREATMENT, DRUG INJECTION (OPTIONAL):

- Rand suggests injecting medication or active molecules through the Clave connectors (blue) integrated at different points in the circuit.
- In order to obtain the best therapeutic result, it is advisable to inject drug only when the perfused district has already reached the required treatment temperature.

WARNING

Every decision about the temperature level, duration of the perfusion and catheter positioning is remitted to the physician responsible of the therapy, as well as the dosage and drug injection method, consistent with the established usage limits of Rand devices.

The compatibility of use of the quoted devices with some commonly used medications - at the discretion of the physician - in intraperitoneal and/or intrapleural hyperthermal perfusion treatments has been appropriately verified by careful review of scientific literature and analysis of clinical test data relating to use of commercial medications based on the following active ingredients of consolidated use: CARBOPLATIN, CISPLATIN, MITOMYCIN C, ADRIAMYCIN/DOXORUBICIN, IRINOTECAN, MELPHALAN and OXALIPLATIN. Rand does not assume any responsibility for use of the devices in combination with medications or preparations not previously evaluated.

END OF TREATMENT LIQUID RINSING AND RECOVERY

Intra-peritoneal and/or intra-pleural hyperthermal perfusion treatments include a recovery and a rinsing phase. At the end of the therapy phase a simultaneous recovery and rinsing phase is started for the removal of the solution used for hyperthermal perfusion and of possible waste products. These steps are accomplished by perfusing the anatomical district with a predetermined volume of fresh washing solution. Due to potential high residual toxicity of the perfusion solutions and to allow the medical staff to perform safe perfusion and rinsing liquids recovery and disposal procedures, follow the indications below:

RINSING PHASE:

- **Before** activating Rinsing phase, replace the solution bags hanging on the IV pole with the new sterile solution bags that will be used for the rinsing procedure.
- **During** Rinsing phase (in order to avoid any balance error during Circulation phase) it is possible to connect the additional 5 litres and 10 liters waste bags (codes R2500155 and R9900111 respectively, included in the kit) to the reservoir bag, in case the volume to recover from patient's body cavity is greater than 8 liters (maximum capacity of the reservoir bag). Connect first the 5 liters waste bag and then the 10 liters waste bag (included in the kit) to the female luer connector at the bottom of the reservoir and open the clamp on the luer connector.

COMPONENTS REPLACEMENT

As a precaution a spare device should always be available during priming and therapy phases. If you believe that a particular event (such as insufficient performances, leaks, filter saturation, etc.) may affect the patient/user safety, interrupt the treatment as explained in the PERFORMER HT equipment operating manual (Section: "Forced Interruption of the Treatment") and proceed to the line replacement.

MEDICAL DEVICES TO BE USED IN COMBINATION

The HANG&GO BASIC HV is specifically designed for use in combination with the RAND PERFORMER HT equipment.

The terminals of the table pack lines (patient inflow and outflow) are compatible with Ch28 catheters/drainages.

The Flow Reverse (code R9900101), if used (optional), is an accessory intended to be connected between the Hang&Go equipment lines and the table pack lines. For set-up and use of this accessory, please refer to the instructions for use provided in the Flow Reverse package.

WARNING

The compatibility of use of the devices listed in this instruction leaflet has been appropriately proved and verified through in vitro and in vivo tests. Rand will not be liable for any consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to Rand or to its distributor and the competent authority of the Member State in which the user is established.

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RAND warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RAND warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RAND cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RAND, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RAND undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RAND up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RAND or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RAND disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RAND, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.