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DESCRIPTION

Disposable device manufactured/assembled by RAND. Not made with natural rubber latex. Made of non-pyrogenic, non-toxic and biocompatible materials according to ISO 10993 standards and in relation with the intended use.

Specifically designed for use in regional cancer therapies for hyperthermic isolated limb perfusion.

TECHNICAL SPECIFICATIONS

Recommended range of blood flow: above 200 ml/min (as displayed on the screen).

The table below shows the correlation between the set flow rate (as seen on the Performer display) and the real flow circulating into the device:

Set flow [ml/min]	Real flow [ml/min]
100	15
200	30
300	45
400	60
500	75
600	90
700	105
800	120
900	135
1000	150
1100	165
1200	180

INTENDED USE

The device is an accessory for the ILP HT line set and it allows to conduct ILP treatment with low flow rates in order to treat small arm portions (e.g. forearm).

SIDE EFFECTS AND CONTRAINDICATIONS

There is no evidence of side effects specifically related to the use of the product. The device is contraindicated for procedures outside of the intended use.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer equipment user's manual.

SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

WARNING

Indicates serious consequences and potential safety hazards for the user and/or patient arising from the normal or improper use of the device, together with the limitations on use and the measures to be taken if such events occur.

CAUTION

Indicates every possible precaution to be taken by the user for the safe and effective use of the device.

Below follows some general safety information for proper use of the device. Furthermore, specific safety-related information is given in the instructions for use where they may condition the operation to be performed.

For further information and/or in case of complaints, contact Rand directly or the authorised local representative.

WARNING

- Carefully read all the information before use. Failure to understand or follow the instructions or reading them without due care may result in serious injury to or the death of the patient.
- The device is intended for professional use. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.
- Carry out a visual inspection and carefully check the devices before use. Transport and/or storage conditions other than those prescribed may have caused damage to the devices.
- Sterility is guaranteed provided the sterile packaging is not damp, open, tampered with or damaged. Do not use the devices if sterility is not guaranteed. Check the use-by date on the label. Do not use after this date.
- The devices must be used immediately after opening the sterile packaging.
- Introduce and maintain the correct dosage of anticoagulant (where required) and monitor this before, during and after treatment.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo further processing. Resterilisation may be ineffective and might compromise the integrity of the device, causing serious harm to the patient.
- Connect the lines in such a way as to prevent the tubing becoming bent or squeezed.
- Avoid alcohol, anaesthetic fluids (i.e. isoflurane), or corrosive solvents (i.e. acetone) coming into contact with the device as they may compromise its structural integrity.
- Because of thermal dissipation, using blood flow rates lower than those recommended does not guarantee that the hyperthermia is attained.

CAUTION

- After use, dispose of the device in compliance with the applicable regulations in force in the country of use.

SETUP

WARNING

This procedure is to be carried out using an aseptic method.

WARNING

The following replacement procedures must be carried out before setting up the ILP HT line set on the equipment.

REPLACING THE PLASTIC HOLDERS OF PM2 PUMP


1. Open PM2 pump cover.
2. Unscrew the screw with black knob under the PM2 pump.
3. Remove the lower front plate with the two white plastic holders for tubing.
4. Replace the holders with those for 4x8 tubing provided by RanD, taking care to position the openings the same way as the previous ones.
5. Put the lower front plate with the new holders back in its housing into PM2 and screw it again.

REPLACING THE PM2 PUMP SEGMENT

CAUTION


It is recommended to execute the cutting and the related connection described in 1 and 2 below on one side of the segment first, and then on the other side, to minimize contamination risks.

1. Take the ILP HT line set R9900086 and cut, with a sterile tool, the pump segment labeled PM2-IN of the OXY INLET LINE R2500184 as indicated in picture 1, i.e. between the pump segment and the connection to the priming line with the YELLOW plastic clamp, and between the pump segment and the connection to the emptying line with the GREEN plastic clamp.

Picture 1

2. Replace the pump segment cut off in point 1 with that of the ILP HT LF accessory, connecting it to the cut points indicated by the arrows in picture 2.

Picture 2

MEDICAL DEVICES THAT CAN BE USED IN COMBINATION

The accessory ILP HT LF has been specifically designed for use in combination with the ILP HT line set code R9900086.

WARNING

The compatibility of the devices listed in the present instructions has been proved and verified with significant in vitro and in vivo tests. Rand will not be liable for any consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to Rand or to its distributor and the competent authority of the Member State in which the user is established

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RanD warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RanD warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RanD cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RanD, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RanD undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RanD up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RanD or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RanD disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RanD, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.