*Instructions for use
Istruzioni per l'uso
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Οδηγίες χρήσης
Gebruiksaanwijzing
Инструкция по использованию
Instrukcja użycia
Uputstvo za upotrebu
Указания за употреба*

*Upute za upotrebu
Пайдалану бойынша нұсқаулар
Instruções de utilização
Instrucţiuni de utilizare
Käyttöohjeet*

R9900086 | ILP HT

EN	<i>Lines set for Limb or Organ Perfusion</i>
IT	<i>Set di linee per Perfusioni d'arto o d'organo</i>
FR	<i>Set de lignes pour Perfusion de Membres ou d'Organes</i>
DE	<i>Leitungsset für Perfusion von Gliedmaßen oder Organen</i>
ES	<i>Juego de líneas para Perfusion de Extremidades u Órganos</i>
EL	<i>Σετ γραμμής για έκπλυση οργάνου ή άκρου</i>
NL	<i>Lijnenset voor de perfusie van ledematen of organen</i>
RU	<i>Система для изолированной гипертермической перфузии конечности или органа</i>
PL	<i>Zestaw linii do perfuzji narządu lub kończyny</i>
SR	<i>Komplet linija za perfuziju uda ili organa</i>
BG	<i>Набор линии за перфузия на крайници или органи</i>
HR	<i>Komplet linija za perfuziju ekstremiteta ili organa</i>
KK	<i>Мүше немесе орган перфузиясы үшін орнатылған желілер</i>
PT	<i>Conjunto de linhas para Perfusão de Membros ou Órgãos</i>
RO	<i>Set de linii pentru perfuzie la nivelul membrelor sau organelor</i>
FI	<i>Letkusarja raajojen tai elinten perfuusiota varten</i>

EN. Do not reuse
IT. Non riutilizzare
FR. Ne pas réutiliser
DE. Nicht zur Wiederverwendung
ES. No reutilizar
TR. Tek kullanımlık
EL. Μίας χρήσης
NL. Niet opnieuw gebruiken
RU. Запрет на повторное применение
PL. Do jednorazowego użytku
SR. Nije za višekratnu upotrebu
LT. Nenaudokite pakartotinai
SV. Får inte återanvändas
BG. Да не се използва повторно
DA. Må ikke genbruges
HR. Nije za višekratnu upotrebu
ZH. 请不要重复使用
ZH. 不可重複使用
CS. Nepoužívejte opakovaně
KK. Қайта пайдалануға болмайды
PT. Não reutilizar
RO. A nu se reutiliza
FI. Ei saa käyttää uudelleen

EN. Content
IT. Contenuto
FR. Contenu
DE. Inhalt
ES. Contenido
TR. İçerik
EL. Περιεχόμενο
NL. Inhoud
RU. Содержит
PL. Zawiera
SR. Sadržaj
LT. Sudėtis
SV. Innehåll
BG. Съдържание
DA. Indhold
HR. Sadržaj
ZH. 内容
ZH. 潤德
CS. Obsah
KK. Содержимое
PT. Conteúdo
RO. Conținut
FI. Sisältö

EN. Non-pyrogenic
IT. Apirogeno
FR. Apyrogène
DE. Pyrogenfrei
ES. Apirogeno
TR. Apirojen
EL. Μη πυρετογόνο
NL. Niet-pyrogeen
RU. Апирогенно
PL. Niepirogenny
SR. Ne stvara toplotu
LT. Nepirogeniškas
SV. Pyrogenfri
BG. Апирогенно
DA. Ikke-pyrogen
HR. Nepirogeno
ZH. 无热源的
ZH. 無致熱原
CS. Nepyrognní
KK. Апирогенді
PT. Não pirogénico
RO. Apirogen
FI. Ei-pyrogeninen

EN. Quantity
IT. Quantità
FR. Quantité
DE. Menge
ES. Cantidad
TR. Adet
EL. Ποσότητα
NL. Aantal
RU. Количество
PL. Ilość
SR. Količina
LT. Kiekis
SV. Antal
BG. Количество
DA. Antal
HR. Količina
ZH. 数量
ZH. 數量
CS. Množství
KK. Саны

EN. Keep dry
IT. Teme l'umidità
FR. Craint l'humidité
DE. Trocken lagern
ES. Matener seco
TR. Kuru saklayın
EL. Διατηρήστε στεγνό
NL. Droog bewaren
RU. Не допускать воздействия влаги
PL. Przechowywać w suchym miejscu
SR. Čuvati na suvom
LT. Laikykite sausai
SV. Förvaras torrt
BG. Да се държи на сухо място
DA. Skal holdes tørt
HR. Čuvati na suhom
ZH. 请保持干燥
ZH. 保持乾燥
CS. Uchovávejte v suchu
KK. Құрғақ сақтаңыз
PT. Manter seco
RO. A se păstra ferit de umezeală
FI. Säilytettävä kuivana

EN. Date of manufacture
IT. Data di fabbricazione
FR. Date de fabrication
DE. Herstellungsdatum
ES. Fecha de fabricación
TR. Üretim tarihi
EL. Ημερομηνία κατασκευής
NL. Fabricagedatum
RU. Дата производства
PL. Data produkcji
SR. Datum proizvodnje
LT. Pagaminimo data
SV. Tillverkningsdatum
BG. Дата на производство
DA. Fremstillingsdato
HR. Datum proizvodnje
ZH. 生产日期
ZH. 製造日期
CS. Datum výroby
KK. Өндүрү күні
PT. Data de fabrico
RO. Data fabricației
FI. Valmistuspäivämäärä

EN. Do not re-sterilize
IT. Non ristilizzare
FR. Ne pas restérilizer
DE. Nicht reesterilisieren
ES. No volver a esterilizar
TR. Tekrar steril etmeyin
EL. Μην αποστειρώνετε ξανά
NL. Niet opnieuw steriliseren
RU. Не стерилизовать повторно
PL. Nie sterylizować
SR. Ne sterilizovati više puta
LT. Nesterilizuoti pakartotinai
SV. Får inte återsteriliseras
BG. Да не се стерилизира повторно
DA. Må ikke gensteriliseres
HR. Ponovna sterilizacija nije dopuštena
ZH. 请勿重复消毒
ZH. 不可重複滅菌
CS. Nesterilizujte opakovaně
KK. Қайта стерилдеуге болмайды
PT. Não reesterilizar
RO. A nu se reesteriliza
FI. Ei saa steriloida uudelleen

EN. Use by
IT. Utilizzare entro
FR. Utiliser jusque
DE. Verwendbar bis
ES. Fecha de caducidad
TR. Son kullanım tarihi
EL. Χρήση έως
NL. Te gebruiken vóór
RU. Использовать до
PL. Zużyć do
SR. Rok upotrebe
LT. Tinka naudoti iki
SV. Använd före
BG. Годен до
DA. Anvendes inden
HR. Upotrijebiti do
ZH. 使用者
ZH. 保存期限
CS. Spotřebujte do
KK. Жарамдылық мерзімі

EN. Batch code
IT. Numero di lotto
FR. Numéro de lot
DE. Chargenbezeichnung
ES. Número de lote
TR. Lot numarası
EL. Κωδικός παρτίδας
NL. Partijnummer
RU. Номер партии
PL. Kod partii
SR. Šifra serije
LT. Partijos numeris
SV. Batchnummer
BG. Партиден номер
DA. Partikode
HR. Šifra serije
ZH. 生产批号
ZH. 製造批號
CS. Kód šarže
KK. Жинақ коды
PT. Código de lote
RO. Număr lot
FI. Eräkoodi

EN. Manufacturer
IT. Fabbricante
FR. Fabricant
DE. Hersteller
ES. Fabricante
TR. Üretici
EL. Κατασκευαστής
NL. Fabrikant
RU. Производитель
PL. Producent
SR. Proizvođač
LT. Gamintojas
SV. Tillverkare
BG. Производител
DA. Producent
HR. Proizvođač
ZH. 生产厂家
ZH. 製造廠
CS. Výrobce
KK. Өндүрүшү
PT. Fabricante
RO. Producător
FI. Valmistaja

EN. Assembled by
IT. Assemblato da
FR. Assemblé par
DE. Zusammengefügt von
ES. Ensamblado por
TR. Birleştirici
EL. Συναρμολόγηση από
NL. Geassembleerd door
RU. Упаковщик
PL. Złożono
SR. Sastavio
LT. Surinko
SV. Monterad av
BG. Сглобено от
DA. Samlet af
HR. Sastavio
ZH. 组装
CS. Sestavil
KK. Құрастырушы
PT. Montado por
RO. Asamblat de
FI. Koonpanija

REF

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■ ENGLISH ■

DESCRIPTION

Sterile single-use devices designed and intended to be used in combination with RanD Performer HT equipment to enable the entire system to reach the intended purpose, i.e. to provide warmed fluids to the patient in the Isolated Limb / Organ Perfusion (also called ILP) in the oncology field.

Specifically, they are extracorporeal circulation tubing sets realized with biocompatible medical grade plastic materials intended to physically connect the patient to the dedicated equipment in order to allow the circulation of the warmed fluids from the equipment to the patient and back, in a continuous closed loop circulation.

TECHNICAL SPECIFICATIONS

Recommended blood flow: from 100 to 1200 ml/min

Water circuit temperature range: from 28 to 45°C

Blood circuit priming volume: 340 ml

Global priming volume: 700 ml (taking into account 200-250 ml in the oxygenator venous reservoir).

INTENDED USE

The disposable circuit for extracorporeal circulation "ILP HT" code R9900086, is intended for use during hyperthermia-based therapies for isolated limb or organ perfusion ("ILP" treatments). The device may be used safely in combination with the medical devices listed in the paragraph entitled "Medical devices that can be used in combination".

SIDE EFFECTS AND CONTRAINDICATIONS

There is no evidence of side effects specifically related to the use of the product.

The device is contraindicated for procedures outside of the intended use.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer equipment user's manual.

SAFETY INFORMATION

Information intended to attract the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

WARNING

Indicates serious consequences and potential safety hazards for the user and/or patient arising from the normal or improper use of the device, together with the limitations on use and the measures to be taken if such events occur.

CAUTION

Indicates every possible precaution to be taken by the user for the safe and effective use of the device.

Some general safety information is given below in order to warn operators about to use the device. Further specific safety information is also given in the instructions for use where this may affect the operation to be carried out.

WARNING



- Carefully read all the information before use. Failure to understand or follow the instructions or reading them without due care may result in serious injury to or the death of the patient.
- The device is intended for professional use. The procedure must be carefully and constantly controlled by qualified personnel trained in use of the device.
- Carry out a visual inspection and carefully check the devices before use. Transport and/or storage conditions other than those prescribed may have caused damage to the devices.
- Sterility is guaranteed provided the sterile packaging is not damp, open, tampered with or damaged. Do not use the devices if sterility is not guaranteed. Check the use-by date on the label. Do not use after this date.
- The devices must be used immediately after opening the sterile packaging.
- Introduce and maintain the correct dosage of anticoagulant (where required) and monitor this before, during and after treatment.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo further processing. Resterilisation may be ineffective and might compromise the integrity of the device, causing serious harm to the patient.
- Connect the lines in such a way as to prevent the tubing becoming bent or squeezed.
- Avoid alcohol, anaesthetic fluids (i.e. isoflurane), or corrosive solvents (i.e. acetone) coming into contact with the device as they may compromise its structural integrity.
- After use, dispose of the device in compliance with the applicable regulations in force in the country of use.
- For further information and/or in the event of a complaint, contact Rand or the authorised local representative.

SETUP

WARNING

This procedure is to be carried out using an aseptic method.

1. Extract the "ILP HT" line set, code R9900086, from its external packaging and remove the sterile pack.
2. Carefully place the pump segments PM1 and PM2 paying attention to the fitting direction; bear in mind that all the pumps in the Performer HT equipment turn clockwise. In order to facilitate the correct set up of the lines, each pump segment is marked with a yellow adhesive collar on the inlet side, which indicates the installation side and the corresponding pump.
3. Other tubing, or tubing portions, to be placed in the equipment automatic clamps are identified by a label marked CL1 or CL2 (according to the fluid diagram required by the therapy) and by their relevant colour code (red or blue).
4. Check the position and integrity of the membrane transducers used on the pressure monitoring lines; Depending on the pressure to be monitored (positive or negative), the deformable membrane may be positioned on one side or the other of the transducer. To facilitate the membrane correct positioning, the body of each transducer bears a yellow label on the side on which the membrane must be positioned for correct use. In addition, the label shows the identification mark of the machine pressure sensor to which the pressure transducer must be connected (from PR1 to PR6). The diagram below shows the correct positioning of the deformable membrane used to monitor:

If necessary, connect a syringe to the fluid isolator membrane connector and then inject or remove air until the membrane is correctly oriented.

5. Connect the pressure lines to the corresponding transducer on the apparatus.
6. Position the heating bag of the HE inlet line, code R2500264, in the correct housing in the heater. Following the equipment on-line instructions, join the quick connector to the corresponding connector placed on the 6000 ml reservoir bag code R2500267, previously placed in the proper holders. Complete assembly of the line by connecting the outlet heater line end to the oxygenator heat exchanger inlet port.
7. Complete the heating circuit by connecting the heat exchanger outlet to the 6000 ml bag, using the HE outlet line, code R2500265.
8. It is advisable to feed the oxygenator through the filtered gas line code R2500266 connected to appropriate gas inlet port.
9. Join the R2500263 oxygenator inlet line ends respectively to the cardiomy reservoir outlet and to the oxygenator inlet, following the instructions provided on the PERFORMER equipment display, close the yellow and green clamps placed on the line.
10. Complete the assembly by connecting the terminals of the A/V line, code R2500262, respectively to the arterial outlet of the oxygenator and to the venous drain located on the venous cardiomy reservoir; if necessary, secure the connections with appropriate plastic locking clips.
11. Hang the bags containing the priming solution on the stand.
12. Insert the spikes of the two lines in the solution bags.
13. Check the correct configuration of the complete circuit.

FILLING PROCEDURE

Follow the step-by-step instructions listed in the operating manual of the Performer HT equipment for a correct handling and set up of the circuit.

CAUTION

- 1) Always use at least 2000 ml of normal saline for the filling procedure.
- 2) During this phase, it is advisable to tap the outside of the devices gently, if necessary, in order to aid the complete evacuation of the air present in the circuit, particularly in the low-flow segments.
- 3) In order to evacuate the air possibly trapped into the oxygenator venous reservoir, it is recommended to open the oxygenator purge/recirculation line (in some models these two lines coincide, in others, the purge line is located at the top of the oxygenator unit and the recirculating line is lower down, near the venous reservoir) for a few seconds, immediately after the end of the priming phase.

STARTING THERAPY

If all the conditions have been met and the instructions followed for completing the priming/rinsing phase correctly, the operator can initiate the therapy by pressing START directly on the equipment touch screen.

CAUTION

- 1) Before starting, check once more that the devices connected to the circuit and/or used in combination with it are correctly filled.

■ ENGLISH ■

- 2) For an adequate therapy monitoring connect the temperature probes cabling to the appropriate jacks on the HTS module of the PERFORMER HT system (from T1 to T8):
 - By using the appropriate cable, connect the temperature probes located on the A/V line to channels T1 and T3 on the HTS module;
 - By using the appropriate cable, connect the temperature probe located on the A/V line to channel T2 of the HTS module;
 - Connect any other temperature probe to the corresponding channels on the HTS module, if applicable;
- 3) Open the external package, release the R2500262 A/V line table pack and clamp both the venous cannulae return line and withdrawal cannulae line. Remove the "U" silicone connector which joins the withdrawal and infusion segment line and connect them respectively to the venous and arterial cannulae previously positioned in the perfusion district vascular accesses.
- 4) When the blood begins to flow into the circuit, start the oxygen flow towards the oxygenator gas inlet. To avoid gas emboli in the blood circuit, make sure that the pressure inside the blood compartment is always higher than the pressure in the gas compartment.
- 5) RAND advises to reach the right treatment parameters gradually, to ensure a correct therapy start.

WARNING

The ILP devices have been designed for use up to a maximum pressure of 500 mmHg. Rand accepts no responsibility whatsoever for improper setting of the parameters or accidental pressurisation of the devices to levels incompatible with the stated limit.

THERAPY

GENERAL INSTRUCTIONS:

For a correct execution of this phase, follow the indications provided with the Performer HT equipment (user's manual and equipment display).

- During the extracorporeal circulation phase it is recommended to continuously monitor the level of the blood drained from the limb / organ venous oxygenator reservoir, and properly adjust the blood flow to avoid the total depletion of the venous cardiotomy reservoir and the consequent possible access of air in the blood circuit.
- Few minutes after the therapy start it is recommended to check the blood gas composition (by means of a hemogasanalyzer); if necessary correct the oxygen concentration of the ventilation gas (FiO2) by turning the mixing flowmeter adjustment knob.
- Constantly monitor the correct dose and efficacy of the adopted anticoagulant treatment, before, during and after the therapy.

DRUG INJECTION (OPTIONAL)

- Rand recommends to inject the drug directly into the oxygenator venous reservoir or through the appropriate clamp-connector fitted to the oxygenator inlet line, code R2500263.
- In order to obtain the best therapeutic result it is recommended to inject the drug only when the perfused limb/organ has already reached the therapy established temperature.

WARNING

Every decision about the temperature level, perfusion duration and cannulae positioning is remitted to the physician responsible of the therapy, as well as the dosage and drug injection method, in accordance with the established usage limits of Rand devices.

The compatibility of use of the quoted devices with some commonly used medications - at the discretion of the physician - in isolated hyperthermic perfusion treatments has been appropriately verified by careful review of scientific literature and analysis of clinical test data relating to use of commercial compounds based on the following active ingredients of consolidated use: TNF- α , ADRIAMYCIN/DOXORUBICIN, OXALIPLATIN, and MELPHALAN. Rand does not assume any responsibility for use of the devices in combination with medications or preparations not previously evaluated.

THERAPY COMPLETION; RECOVERY OF FLUIDS

"LRT" therapies contemplate the recovery of the fluids used during the hyperthermic perfusion, which may be used in combination with drugs. Due to their high residual toxicity level, the perfusion fluids must be physically isolated and collected separately. To this purpose the "ILP" circuits provide a special system for the recovery and management of the fluids used during the therapy. The entire volume of such fluids, in the worst-case scenario, can be collected in suitable waste bags (code R2500268) that allow a safe management of the emptying phase.

RINSING//REINFUSION PHASE:

Before activating this phase:

- Hang the solution bag to be used for the rinsing procedure on the stand;
- Connect the additional 5 litre waste bag (code R2500268, included in the set) to the connector fitted to the A/V line, code R2500262.

PARTS REPLACEMENT

As a precaution, a spare device must always be available during the priming phase and in the course of the therapy.

■ ENGLISH ■

If you feel that a particular event (such as insufficient performances, leaks, abnormal blood parameters, et c.) might affect the patient/user safety, proceed to the line replacement as explained in the PERFORMER HT equipment operating manual.

MEDICAL DEVICES THAT CAN BE USED IN COMBINATION

The ILP HT line set, code R9900086, has been specifically designed for use in combination with the RAND PERFORMER HT apparatus.

The line set has been optimised for connection to oxygenators whose specifications are compatible with those set out in the paragraph "technical specifications" and which are fitted with the following connectors: 1/4" connector for blood inlet and outlet, 3/8" connector for rigid venous reservoir and Hansen-type connector for the heat exchanger.

WARNING

The compatibility of the devices listed in the present instructions has been proved and verified with significant in vitro and in vivo tests. Rand will not be liable for any consequences, problems or malfunctions deriving from the use of the system in conjunction with different devices not previously evaluated.

CAUTION

The diameter of the circuit connections must be compatible with the dimensions of the connectors assembled on the devices used in combination.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to Rand or to its distributor and the competent authority of the Member State in which the user is established

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RAND warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RAND warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RAND cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RAND, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RAND undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RAND up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RAND or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RAND disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RAND, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.