*Instructions for use
Istruzioni per l'uso
Mode d'emploi
Instrucciones de uso
Brugsanvisninger
Bruksanvisning
Instrukcja użycia
Указания за употреба
Upute za upotrebu*

*Пайдалану бойынша нұсқаулар
Instruções de utilização
Instrucțiuni de utilizare
تعليمات الاستخدام
Használati utasítás
Käyttöohjeet
Petunjuk penggunaan*

R9900089

HYPER PROTECT-OR

EN	<i>Abdominal cover</i>
IT	<i>Copertura addominale</i>
FR	<i>Couverture abdominale</i>
ES	<i>Cobertura abdominal</i>
DA	<i>Abdominal dækken</i>
SV	<i>Bukhölje</i>
PL	<i>Ośłona na jamę brzuszną</i>
BG	<i>Абдоминална покривка</i>
HR	<i>Abdominalni pokrov</i>
KK	<i>Құрсақ жабыны</i>
PT	<i>Tampa abdominal</i>
RO	<i>Protecție abdominală</i>
AR	<i>غطاء البطن</i>
HU	<i>Hasi zárótakaró</i>
FI	<i>Vatsansuojus</i>
ID	<i>Penutup perut</i>

EN. Do not reuse
IT. Non riutilizzare
FR. Ne pas réutiliser
DE. Nicht zur Wiederverwendung
ES. No reutilizar
TR. Tek kullanımlık
NL. Niet opnieuw gebruiken
RU. Не использовать повторно
PL. Do jednorazowego użytku
SR. Nije za višekratnu upotrebu
LT. Nenaudokite pakartotinai
SV. Får inte återanvändas
BG. Да не се използва повторно
DA. Må ikke genbruges
HR. Nije za višekratnu upotrebu
ZH. 请不要重复使用
ZH. 不可重复使用
CS. Nepoužívejte opakovaně
KK. Қайта пайдалануға болмайды
PT. Não reutilizar
RO. A nu se reutiliza
AR. لا تُعد استخداماً
HU. Ne használja újra
FI. Ei saa käyttää uudelleen
ID. Jangan gunakan kembali

EN. Content
IT. Contenuto
FR. Contenu
DE. Inhalt
ES. Contenido
TR. İçerik
NL. Inhoud
RU. Содержит
PL. Zawiera
SR. Sadržaj
LT. Sudėtis
SV. Innehåll
BG. Съдържание
DA. Indhold
HR. Sadržaj
ZH. 内容
ZH. 潤德
CS. Obsah
KK. Содержимое
PT. Conteúdo
RO. Conținut
AR. المحتوى
HU. Tartalom
FI. Sisältö
ID. Isi

EN. Non-pyrogenic
IT. Apirogeno
FR. Apyrogène
DE. Pyrogenfrei
ES. Apirogeno
TR. Apirojen
NL. Niet-pyrogeen
RU. Апирогенно
PL. Niepirogenny
SR. Ne stvara toplotu
LT. Nepirogeniškas
SV. Pyrogenfri
BG. Апирогенно
DA. Ikke-pyrogen
HR. Nepirogeno
ZH. 无热源的
ZH. 無致熱原
CS. Nepyrogní
KK. Апирогенді
PT. Não pirogénico
RO. Apirogen
AR. غير بيريوجيني
HU. Nem pirogén
FI. Ei-pyrogeeninen
ID. Non-pirogenik

EN. Quantity
IT. Quantità
FR. Quantité
DE. Menge
ES. Cantidad
TR. Adet
NL. Aantal
RU. Количество
PL. Ilość
SR. Količina
LT. Kiekis
SV. Antal
BG. Количество
DA. Antal
HR. Količina

EN. Keep dry
IT. Teme l'umidità
FR. Craint l'humidité
DE. Trocken lagern
ES. Matener seco
TR. Kuru saklayın
NL. Droog bewaren
RU. Держать сухим
PL. Przechowywać w suchym miejscu
SR. Čuvati na suvom
LT. Laikykite sausai
SV. Förvaras torrt
BG. Да се държи на сухо място
DA. Skal holdes tørt
HR. Čuvati na suhom
ZH. 请保持干燥
ZH. 保持乾燥
CS. Uchovávejte v suchu
KK. Құрақ сақтаңыз
PT. Manter seco
RO. A se păstra ferit de umezeală
AR. حافظ على جفافه
HU. Tartsa szárazon
FI. Säilytettävä kuivana
ID. Jaga agar tetap kering

EN. Date of manufacture
IT. Data di fabbricazione
FR. Date de fabrication
DE. Herstellungsdatum
ES. Fecha de fabricación
TR. Üretim tarihi
NL. Fabricagedatum
RU. Дата производства
PL. Data produkcji
SR. Datum proizvodnje
LT. Pagaminimo data
SV. Tillverkningsdatum
BG. Дата на производство
DA. Fremstillingsdato
HR. Datum proizvodnje
ZH. 生产日期
ZH. 製造日期
CS. Datum výroby
KK. Өндіру күні
PT. Data de fabrico
RO. Data fabricației
AR. تاريخ التصنيع
HU. Gyártás dátuma
FI. Valmistuspäivämäärä
ID. Tanggal pembuatan

EN. Do not re-sterilize
IT. Non ristilizzare
FR. Ne pas restériliser
DE. Nicht reesterilisieren
ES. No volver a esterilizar
TR. Tekrar steril etmeyin
NL. Niet opnieuw steriliseren
RU. Не стерилизовать повторно
PL. Nie sterylizować
SR. Ne sterilizovati više puta
LT. Nesterilizuokite pakartotinai
SV. Får inte återsteriliseras
BG. Да не се стерилизира повторно
DA. Må ikke gensteriliseres
HR. Ponovna sterilizacija nije dopuštena
ZH. 请勿重复消毒
ZH. 不可重复滅菌
CS. Nesterilizujte opakovaně
KK. Қайта стерилдеуге болмайды
PT. Não reesterilizar
RO. A nu se reesteriliza
AR. لا تُعد تعقيم
HU. Ne sterilizálja újra
FI. Ei saa steriloida uudelleen
ID. Jangan disterilkan ulang

EN. Use by
IT. Utilizzare entro
FR. Utiliser jusqu'à
DE. Verwendbar bis
ES. Fecha de caducidad
TR. Son kullanım tarihi
NL. Te gebruiken vóór
RU. Использовать до
PL. Zużyć do
SR. Rok upotrebe
LT. Tinka naudoti iki
SV. Använd före
BG. Годен до
DA. Anvendes inden
HR. Upotrijebiti do

EN. Batch code
IT. Numero di lotto
FR. Numéro de lot
DE. Chargenbezeichnung
ES. Número de lote
TR. Lot numarası
NL. Partijnummer
RU. Номер партии
PL. Kod partii
SR. Šifra serije
LT. Partijos numeris
SV. Batchnummer
BG. Партиден номер
DA. Partikode
HR. Šifra serije
ZH. 生产批号
ZH. 製造批號
CS. Kód šarže
KK. Жинақ коды
PT. Código de lote
RO. Număr lot
AR. رمز التشغيل
HU. Tételkód
FI. Eräkoodi
ID. Kode batch

EN. Manufacturer
IT. Fabbricante
FR. Fabricant
DE. Hersteller
ES. Fabricante
TR. Üretici
NL. Fabrikant
RU. Производитель
PL. Producent
SR. Proizvođač
LT. Gamintojas
SV. Tillverkare
BG. Производител
DA. Producent
HR. Proizvođač
ZH. 生产厂家
ZH. 製造廠
CS. Výrobce
KK. Өндіруші
PT. Fabricante
RO. Producător
AR. الشركة المصنعة
HU. Gyártó
FI. Valmistaja
ID. Produsen

EN. Assembled by
IT. Assemblato da
FR. Assemblé par
DE. Zusammengefügt von
ES. Ensamblado por
TR. Birleştirici
NL. Geassembleerd door
RU. Упаковщик
PL. Złożono
SR. Sastavio
LT. Surinko
SV. Monterad av
BG. Сглобено от
DA. Samlet af
HR. Sastavio
ZH. 组装
CS. Sestavil
KK. Құрастырушы
PT. Montado por
RO. Asamblat de
AR. التجميع من قبل
HU. Összeállította
FI. Kokoonpanija
ID. Dirakit oleh

ZH. 数量
ZH. 數量
CS. Množství
KK. Саны
PT. Quantidade
RO. Cantitate
AR. الكمية
HU. Mennyiség
FI. Määrä
ID. Kuantitas

EN. Catalogue number
IT. Numero di catalogo
FR. Référence du catalogue
DE. Bestellnummer
ES. Número de catálogo
TR. Referans numarası
NL. Catalogusnummer
RU. Номер по каталогу
PL. Numer katalogowy
SR. Kataloški broj
LT. Katalogo numeris
SV. Katalognummer
BG. Каталожен номер
DA. Katalognummer
HR. Broj u katalogu
ZH. 目录号
ZH. 產品型號
CS. Katalogové číslo
KK. Каталог нөмірі
PT. Número de catálogo
RO. Număr de catalog
AR. رقم الكتالوج
HU. Katalógusszám
FI. Luettelonumero
ID. Nomor katalog

EN. Keep away from heat
IT. Tenere lontano da sorgenti di calore
FR. Craint la chaleur
DE. Vor Wärmequellen fernhalten
ES. Mantener lejos de fuentes de calor
TR. Isı kaynaklarından uzak tutun
NL. Uit de buurt van warmtebronnen houden
RU. Беречь от тепла
PL. Przechowywać z dala od źródeł ciepła
SR. Držati dalje od toplote
LT. Saugokite nuo šilumos
SV. Skydda från värme
BG. Да се пази от топлина
DA. Må ikke udsættes for varme
HR. Držati dalje od topline
ZH. 远离高温
ZH. 避免日照
CS. Chraňte před teplem
KK. Ыстыктан алыс ұстаңыз
PT. Manter afastado do calor
RO. A se păstra departe de sursele de căldură
AR. أبقيه بعيداً عن الحرارة
HU. Tartsa távol hőtől
FI. Suojattava lämmöltä
ID. Jauhkan dari panas

STERILE EO

EN. Sterile - Ethylene oxide sterilized
IT. Sterile - Sterilizzato ad ossido d'etilene
FR. Stérile - Stérilisé à l'oxyde d'éthylène
DE. Steril - Mit Äthylenoxid sterilisiert
ES. Estéril - Esterilizado por óxido de etileno
TR. Steril - Etilen Oksit
NL. Steriel - Gesteriliseerd met ethyleenoxide
RU. Стерильно - Стерилизовано окисью этилена
PL. Sterylne - Wysterylizowane tlenkiem etylenu
SR. Sterilno - Sterilisano etilen oksidom
LT. Sterilu - Sterilizuota etileno oksidu
SV. Steril - Steriliserad med etylenoxid
BG. Стерилно - Стерилизирано с этиленов оксид
DA. Steril - Steriliseret med ethylenoxid
HR. Sterilno - sterilizirano etilen-oksidiom
ZH. 消毒 - 环氧乙烷消毒法
ZH. EO滅菌
CS. Sterilní - Sterilizováno ethylenoxidem
KK. Стерильді - Этилен тотығымен стерильденген
PT. Estéril - Esterilizado por óxido de etileno
RO. Steril - Sterilizat cu oxid de etilenă
AR. معقم - معقم بأكسيد الإيثيلين
HU. Steril - Etilén-oxiddal sterilizált
FI. Steriili - steriloitu etyleenioksidilla
ID. Steril - Etilen oksida disterilkan

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DESCRIPTION

Hyper Protect-OR is a device intended to delimitate the surgical field surrounding the laparotomy during Hyperthermic Intraperitoneal Perfusion treatments with open/semi-closed abdomen technique. The device reduces heat dispersion and it facilitates reaching and maintaining the temperature in the peritoneal cavity. The device must be placed on a retractor with dimensions compatible with those of the device conical portion (see below).

The corrugated tube available in the package can be connected to the swivel connector to evacuate water vapours and aerosols from the abdomen. The 22F – 15M adapter provided in the package allows the connection of the corrugated tube to aspiration systems compatible with the 22F size of the adapter. An opening placed at the top of the conical portion allows the surgeon's arm to access the cavity.

TECHNICAL FEATURES

Dimensions: 100 x 80 cm.

Dimensions of the conical portion: 60 x 40 x 20 cm.

Wall thickness: 0.3 mm.

Material: PVC.

Tube length: 3 m.

The device is entirely made of medical grade, atoxic and latex-free materials.

INTENDED USE

Sterile disposable device intended to delimitate the surgical field surrounding the laparotomy during Hyperthermic Intraperitoneal Perfusion treatments with open/semi-closed abdomen technique, in order to preserve the target temperature inside the abdominal cavity.

SIDE EFFECTS AND CONTRAINDICATIONS

The device must not be used for procedures outside of the intended use. There is no evidence of side effects specifically related to the use of the product.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer equipment user's manual.

SAFETY INFORMATION

The information to draw the attention of the user to potentially dangerous situations and to ensure correct and safe use of the device are indicated in the text as follows:

WARNING

Indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur both during proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such case.

PRECAUTION

Indicates any special precaution than the user must exercise for safe and effective use of the device.

For further information and/or in case of complaints, contact RanD directly or the authorised local representative.

WARNING

- The user must read and understand all user information prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- The device is intended for professional use only. The procedure must be carefully and constantly controlled by qualified personnel trained in the use of the device.

- The device does not provide a hermetic seal and is not intended as a protective device for personnel. It does not replace the usual personal protective devices recommended during the procedure.
- Visually inspect and carefully check the device before use. Transport and/or storage conditions not conform to the specifications may have caused damage to the devices.
- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the devices if sterility cannot be guaranteed.
- Check the expiry date shown on the label attached. Do not use the devices after this date.
- Use the device immediately after opening the sterile packaging.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo any further processing. Resterilisation may not be effective and may compromise the integrity of the device thus resulting in unintended injury.
- Avoid alcohol, anaesthetic fluids, or corrosive solvents from coming into contact with the devices as they may jeopardize the structural integrity.

PRECAUTION

- After use, dispose of the devices in accordance with the applicable regulations in force in the country where the device is used.

INSTALLATION

Perform the procedure using an aseptic technique.

1. After preparing the Patient for the perfusion (i.e. once placed and connected catheters and temperature probes to the Performer system), remove the Hyper Protect-OR from its package and spread it completely.
2. Orientate the longest side of the abdominal cover of the Hyper Protect-OR along the Patient axis. In order to make hands access to the peritoneal cavity easier, it is advisable to place the cover so that the connection for the tube lays on the opposite side to the operator.
3. Put the conical portion on the Thompson retractor (or equivalent systems with similar features and dimensions) and place it over the surgical field, paying attention to keep it raised from the Patient.
4. Cover the abdomen.
5. In order to reduce possible formation of water vapours and evacuate aerosols, one end of the tube can be connected to the swivel connector on the cover and the other end to a smokes/aerosols aspirator, featuring protective filters with antibacterial and hydrophobic membrane.

MEDICAL DEVICES TO BE USED IN COMBINATION

Hyper Protect-OR can be used with Thompson retractor or with equivalent devices (with similar features and dimensions) and can be connected with most of the vapours and smokes aspiration systems.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to RanD or to its distributor and the competent authority of the Member State in which the user is established

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RAND warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RAND warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RAND cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RAND, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RAND undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RAND up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RAND or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RAND disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RAND, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.