Instructions for use
Istruzioni per l'uso
Mode d'emploi
Gebrauchsanweisung
Instrucciones de uso
Kullanım Talimatları
Οδηγίες χρήσης
Gebruiksaanwijzing
Инструкция по использованию
Instrukcja użytkowania
Uputstvo za upotrebu
Naudojimo instrukcijos
Bruksanvisning

Указания за употреба
Brugsanvisninger
Upute za upotrebu
使用説明書
使用说明
Návod k použití
Паўдалану бойынша нұсқаулар
Instruções de utilização
Instrucţiuni de utilizare
تعليمات الاستخدام
Használati utasítás
Käyttöohjeet
Bruksanvisning
Petunjuk penggunaan

R9900071

EN	Ch24 round silicone Catheter
IT	Catetere tondo in Silicone Ch24
FR	Cathéter arrondi en Silicone Ch24
DE	Runder Silikonkatheter Ch24
ES	Catéter redondo de silicona Ch24
TR	Silikon Ch24 dairesel Kateter
EL	Στρογγυλός καθετήρας από σιλίκονη Ch24
NL	Ronde siliconen katheter Ch24
RU	Круглый силиконовый катетер Ch24
PL	Cewnik okrągły silikonowy Ch24
SR	Okrugli silikonski kateter Ch24
LT	Ch24 apvalus silikoninis kateteris
SV	Ch24 rund silikonkateter
BG	Ch24 кръгъл силиконов катетър
DA	Ch24 rundt silikonekateter
HR	Okrugli silikonski kateter Ch24
ZH	保福腹腔溫熱灌注導管
ZH	Ch24 圓硅胶导管
CS	Okrouhlý silikonový katétr Ch24
KK	Ch24 дөңгелек силикон катетер
PT	Cateter redondo em silicone Ch24
RO	Cateter rotund din silicon Ch24
AR	Ch24 قسطرة سيليكون دائرية
HU	24 Ch-es kerek keresztmetszetű szilikon katéter
FI	Pyöreä koon Ch24 silikonikatetri
NO	Ch24 rundt silikonkateter
ID	Kateter silikon bulat Ch24

EN. Do not reuse
IT. Non riutilizzare
FR. Ne pas réutiliser
DE. Nicht zur Wiederverwendung
ES. No reutilizar
TR. Tek kullanımlık
EL. Μίας χρήσης
NL. Niet opnieuw gebruiken
RU. Не использовать повторно
PL. Do jednorazowego użytku
SR. Nije za višekratnu upotrebu
LT. Nenaudokite pakartotinai
SV. Får inte återanvändas
BG. Да не се използва повторно
DA. Må ikke genbruges
HR. Nije za višekratnu upotrebu
ZH. 请不要重复使用
ZH. 不可重覆使用
CS. Nepoužívejte opakovaně
KK. Кайта пайдалануға болмайды
PT. Não reutilizar
RO. A nu se reutiliza
AR. لا تُعد استخداما
HU. Ne használja újra
FI. Ei saa käyttää uudelleen
NO. Ikke bruk på nytt
ID. Jangan gunakan kembali

EN. Length
IT. Lunghezza
FR. Longueur
DE. Länge
ES. Largo
TR. Uzunluk
EL. Μήκος
NL. Lengte
RU. длина
PL. Długość
SR. Dužina
LT. Ilgis
SV. Längd
BG. Дължина
DA. Længde
HR. Duljina
ZH. 长度
ZH. 長度
CS. Délka
KK. Ұзындығы
PT. Comprimento
RO. Lungime
AR. الطول
HU. Hossz
FI. Pituus
NO. Lengde
ID. Panjang

EN. Non-pyrogenic
IT. Apirogeno
FR. Apyrogène
DE. Pyrogenfrei
ES. Apirogeno
TR. Apirojen
EL. Μη πυρογόνο
NL. Niet-pyrogeen
RU. Апирогенно
PL. Niepirogenny
SR. Ne stvara toplotu
LT. Nepirogeniškas
SV. Pyrogenfri
BG. Апирогенно
DA. Ikke-pyrogen
HR. Nepirogeno
ZH. 无热源的
ZH. 無致熱原
CS. Nepyrognní
KK. Апирогенді
PT. Não pirogénico
RO. Apirogen
AR. غير بيريوجيني
HU. Nem pirogén
FI. Ei-pyrogeenininen
NO. Ikke-pyrogen
ID. Non-pirogenik

EN. Quantity
IT. Quantità
FR. Quantité
DE. Menge
ES. Cantidad
TR. Adet
EL. Ποσότητα
NL. Aantal
RU. Количество
PL. Ilość

EN. Keep dry
IT. Teme l'umidità
FR. Craint l'humidité
DE. Trocken lagern
ES. Matener seco
TR. Kuru saklayın
EL. Διατηρήστε στεγνό
NL. Droog bewaren
RU. Держать сухим
PL. Przechowywać w suchym miejscu
SR. Čuvati na suvom
LT. Laikykite sausai
SV. Förvaras torrt
BG. Да се държи на сухо място
DA. Skal holdes tørt
HR. Čuvati na suhom
ZH. 请保持干燥
ZH. 保持乾燥
CS. Uchovávejte v suchu
KK. Құрғақ сақтаңыз
PT. Manter seco
RO. A se păstra ferit de umezeală
AR. حافظ على جفافه
HU. Tartsa szárazon
FI. Säilytettävä kuivana
NO. Holdes tørt
ID. Jaga agar tetap kering

EN. Date of manufacture
IT. Data di fabbricazione
FR. Date de fabrication
DE. Herstellungsdatum
ES. Fecha de fabricación
TR. Üretim tarihi
EL. Ημερομηνία κατασκευής
NL. Fabricagedatum
RU. Дата производства
PL. Data produkcji
SR. Datum proizvodnje
LT. Pagaminimo data
SV. Tillverkningsdatum
BG. Дата на производство
DA. Fremstillingsdato
HR. Datum proizvodnje
ZH. 生产日期
ZH. 製造日期
CS. Datum výroby
KK. Өндүрүш күні
PT. Data de fabrico
RO. Data fabricației
AR. تاريخ التصنيع
HU. Gyártás dátuma
FI. Valmistuspäivämäärä
NO. Produksjonsdato
ID. Tanggal pembuatan

EN. Do not re-sterilize
IT. Non risterrilizzare
FR. Ne pas restérilizer
DE. Nicht reesterilisieren
ES. No volver a esterilizar
TR. Tekrar steril etmeyin
EL. Μην αποστειρώνετε ξανά
NL. Niet opnieuw steriliseren
RU. Не стерилизовать повторно
PL. Nie sterylizować
SR. Ne sterilizovati više puta
LT. Nesterilizuoti pakartotinai
SV. Får inte återsteriliseras
BG. Да не се стерилизира повторно
DA. Må ikke gensteriliseres
HR. Ponovna sterilizacija nije dopuštena
ZH. 请勿重复消毒
ZH. 不可重覆滅菌
CS. Nesterilizujte opakovaně
KK. Кайта стерилдеуге болмайды
PT. Não reesterilizar
RO. A nu se resteriliza
AR. لا تُعد تعقيم
HU. Ne sterilizálja újra
FI. Ei saa steriloida uudelleen
NO. Ikke steriliser på nytt
ID. Jangan disterilkan ulang

EN. Use by
IT. Utilizzare entro
FR. Utiliser jusque
DE. Verwendbar bis
ES. Fecha de caducidad
TR. Son kullanım tarihi
EL. Χρήση έως
NL. Te gebruiken vóór
RU. Использовать до
PL. Zużyć do

EN. Batch code
IT. Numero di lotto
FR. Numéro de lot
DE. Chargenbezeichnung
ES. Número de lote
TR. Lot numarası
EL. Κωδικός παρτίδας
NL. Partijnummer
RU. Номер партии
PL. Kod partii
SR. Šifra serije
LT. Partijos numeris
SV. Batchnummer
BG. Партиден номер
DA. Partikode
HR. Šifra serije
ZH. 生产批号
ZH. 製造批號
CS. Kód šarže
KK. Жинақ коды
PT. Código de lote
RO. Număr lot
AR. رمز التشغيل
HU. Tételkód
FI. Eräkoodi
NO. Partikode
ID. Kode batch

EN. Manufacturer
IT. Fabbricante
FR. Fabricant
DE. Hersteller
ES. Fabricante
TR. Üretici
EL. Κατασκευαστής
NL. Fabrikant
RU. Производитель
PL. Producent
SR. Proizvođač
LT. Gamintojas
SV. Tillverkare
BG. Производител
DA. Producent
HR. Proizvođač
ZH. 生产厂家
ZH. 製造廠
CS. Výrobce
KK. Өндүрүш
PT. Fabricante
RO. Producător
AR. الشركة المصنعة
HU. Gyártó
FI. Valmistaja
NO. Produsent
ID. Produsen

EN. Size
IT. Dimensioni
FR. Dimensions
DE. Größe
ES. Dimensiones
TR. Boyut
EL. Μέγεθος
NL. Afmetingen
RU. Размер
PL. Rozmiar
SR. Veličina
LT. Dydis
SV. Storlek
BG. Размер
DA. Størrelse
HR. Veličina
ZH. 尺寸
ZH. 规格
CS. Velikost
KK. Өлшемі
PT. Dimensões
RO. Dimensiuni
AR. المقاس
HU. Méret
FI. Koko
NO. Størrelse
ID. Ukuran

EN. Catalogue number
IT. Numero di catalogo
FR. Référence du catalogue
DE. Bestellnummer
ES. Número de catálogo
TR. Referans numarası
EL. Αριθμός καταλόγου
NL. Catalogusnummer
RU. Номер по каталогу
PL. Numer katalogowy

SR.	Količina
LT.	Kiekis
SV.	Antal
BG.	Количество
DA.	Antal
HR.	Količina
ZH.	数量
ZH.	數量
CS.	Množství
KK.	Саны
PT.	Quantidade
RO.	Cantitate
AR.	الكمية
HU.	Mennyiség
FI.	Määrä
NO.	Antall
ID.	Kuantitas

EN.	Not made with natural rubber latex
IT.	Non fatto con lattice di gomma naturale
FR.	Latex exempt
DE.	Latex-frei
ES.	Libre de látex
TR.	Lateks içermez
EL.	Χωρίς λατέξ
NL.	Latexvrij
RU.	Не содержит латекса
PL.	Nie zawiera lateksu
SR.	Ne sadrži lateks
LT.	Be latekso
SV.	Latexfri
BG.	Не съдържа латекс
DA.	Latex-fri
HR.	Bez lateksa
ZH.	请与乳胶隔离
ZH.	非乳膠製品
CS.	Neobsahuje latex
KK.	Латексиз
PT.	Não contém látex
RO.	Fără latex
AR.	خالي من اللاتكس
HU.	Nem természetes gumilatexből készült
FI.	Ei sisällä luonnonkumilateksia
NO.	Ikke laget med naturgummilateks
ID.	Tidak terbuat dari lateks karet alam

EN.	Keep away from heat
IT.	Tenere lontano da sorgenti di calore
FR.	Craint la chaleur
DE.	Vor Wärmequellen fernhalten
ES.	Mantener lejos de fuentes de calor
TR.	Isı kaynaklarından uzak tutun
EL.	Φυλάξτε το μακριά από θερμότητα
NL.	Uit de buurt van warmtebronnen houden
RU.	Беречь от тепла
PL.	Przechowywać z dala od źródeł ciepła
SR.	Držati dalje od toplote
LT.	Saugokite nuo šilumos
SV.	Skydda från värme
BG.	Да се пази от топлина
DA.	Må ikke udsættes for varme
HR.	Držati dalje od topline
ZH.	远离高温
ZH.	避免日照
CS.	Chraňte před teplem
KK.	Ыстықтан алыс ұстаңыз
PT.	Manter afastado do calor
RO.	A se păstra departe de sursele de căldură
AR.	أبقه بعيداً عن الحرارة
HU.	Tartsa távol hőtől
FI.	Suojattava lämmöltä
NO.	Holdes unna varme
ID.	Jauhkan dari panas

STERILE	EO
---------	----

EN.	Sterile - Ethylene oxide sterilized
IT.	Sterile - Sterilizzato ad ossido d'etilene
FR.	Sterile - Stérilisé à l'oxyde d'éthylène
DE.	Steril - Mit Äthylenoxid sterilisiert
ES.	Estéril - Esterilizado por óxido de etileno
TR.	Steril - Etilen Oksit
EL.	Αποστειρωμένο - Αποστειρωμένο με οξείδιο του αιθυλενίου
NL.	Steriel - Gesteriliseerd met ethyleenoxide
RU.	Стерильно - Стерилизовано окисью этилена
PL.	Sterylny - Wysterylizowane tlenkiem etylenu
SR.	Sterilno - Sterilizirano etilen oksidom
LT.	Sterilu - Sterilizuota etileno oksidu
SV.	Steril - Steriliserad med etylenoxid
BG.	Стерилно - Стерилизирано с этиленов оксид
DA.	Steril - Steriliseret med ethylenoxid
HR.	Sterilno - sterilizirano etilen-oksidiom
ZH.	消毒 - 环氧乙烷消毒法

SR.	Rok upotrebe
LT.	Tinka naudoti iki
SV.	Använd före
BG.	Годен до
DA.	Anvendes inden
HR.	Upotrijebiti do
ZH.	使用者
ZH.	保存期限
CS.	Spotřebujte do
KK.	Жарамдылық мерзімі
PT.	Data de validade
RO.	A se utiliza înainte de
AR.	يستخدم من قبل
HU.	Felhasználhatóság dátuma
FI.	Viimeinen käyttöpäivä
NO.	Brukes innen
ID.	Gunakan sampai dengan

EN.	Radioopac
IT.	Radiopaco
FR.	Radio-opaque
DE.	Strahlenundurchlässig
ES.	Radiopaco
TR.	Radyopak
EL.	Ακτινοοσκιερό
NL.	Radiopaak
RU.	Рентгеноконтрастно
PL.	Nieprzepuszczalny dla promieniowania
SR.	Ne propušta zračenje
LT.	Rentgeno spinduliuose kontrastuojantis
SV.	Röntgentät
BG.	Рентгеноконтрастно
DA.	Røntgenfast
HR.	Ne propušta zračenje
ZH.	防辐射
ZH.	輻射無法穿透
CS.	Rentgenkontrastní
KK.	Контрастты зат
PT.	Radiopaco
RO.	Radioopac
AR.	راديوباك
HU.	Sugárfogó
FI.	Röntgenläpinäkyvä
NO.	Radioopac
ID.	Radioopac

EN.	Caution, consult accompanying documents
IT.	Attenzione, vedere le istruzioni per l'uso
FR.	Attention voir notice d'instructions
DE.	Achtung, Begleiddokumente beachten
ES.	Atención, ver instrucciones de uso
TR.	Dikkat, kullanım talimatlarına bakın
EL.	Προσοχή, συμβουλευτείτε τα συνοδευτικά έγγραφα
NL.	Attentie, de bijgevoegde documenten raadplegen
RU.	Внимание! Ознакомиться с сопроводительными документами
PL.	Uwaga! Zapoznać się z dołączonymi dokumentami
SR.	Pažnja, pogledajte prateću dokumentaciju
LT.	Dėmesio, perskaitykite priedamus dokumentus
SV.	Varning – läs de medföljande dokumenten
BG.	Внимание, направете справка в придружаващата документация
DA.	Forsigtig, rådfør dig med de medfølgende dokumenter
HR.	Oprez, pogledati prateće dokumente
ZH.	注意. 请阅读所附文件
ZH.	注意事項，請參考隨附文件
CS.	Upozornění: prostudujte si průvodní dokumentaci
KK.	Ескерту, бірге келген құжаттарды қараңыз
PT.	Cuidado, consultar os documentos fornecidos
RO.	Atenție, consultați documentele însoțitoare
AR.	احترس، راجع الوثائق المرفقة
HU.	Figyelem, olvassa el a mellékelt dokumentumokat
FI.	Huomio, katso liiteasiakirjat
NO.	Forsiktig, se medfølgende dokumenter
ID.	Perhatian, baca dokumen yang disertakan

EN.	Do not use if package is damaged
IT.	Non utilizzare se la confezione è danneggiata
FR.	Ne pas utiliser si l'emballage est endommagé
DE.	Inhalt beschädigter Packung nicht verwenden
ES.	No usar si el paquete está dañado
TR.	Paket hasar görmüşse kullanmayın
EL.	Μην το χρησιμοποιείτε εάν έχει καταστραφεί η συσκευασία
NL.	Niet gebruiken als de verpakking beschadigd is
RU.	Не использовать, если упаковка повреждена
PL.	Nie używać, jeśli opakowanie jest uszkodzone
SR.	Ne koristiti ukoliko je pakovanje oštećeno
LT.	Nenaudokite, jeigu pažeista pakuotė
SV.	Använd inte om förpackningen är skadad
BG.	Да не се използва, ако опаковката е повредена
DA.	Må ikke anvendes, hvis emballagen er beskadiget
HR.	Ne upotrebljavati ako je pakiranje oštećeno
ZH.	包装出现损坏的情况下请不要使用
ZH.	包裝破損不可使用

SR.	Kataloški broj
LT.	Katalogo numeris
SV.	Katalognummer
BG.	Каталоген номер
DA.	Katalognummer
HR.	Broj u katalogu
ZH.	目录号
ZH.	產品型號
CS.	Katalogové číslo
KK.	Каталог нөмірі
PT.	Número de catálogo
RO.	Număr de catalog
AR.	رقم الكتالوج
HU.	Katalógusszám
FI.	Luettelonnumero
NO.	Katalognummer
ID.	Nomor katalog

EN.	Temperature limitation
IT.	Limiti di temperatura
FR.	Limites de température
DE.	Zulässiger Temperaturbereich
ES.	Límite de temperatura
TR.	Sıcaklık limitleri
EL.	Περιορισμός θερμοκρασίας
NL.	Temperatuurlimieten
RU.	Пределы температур
PL.	Ograniczenia temperatury
SR.	Temperaturno ograničenje
LT.	Temperatūros apribojimas
SV.	Temperaturbegränsning
BG.	Температурно ограничение
DA.	Temperaturbegrænsning
HR.	Ograničenje temperature
ZH.	温度限制
ZH.	溫度限制
CS.	Teplotní omezení
KK.	Температура шери
PT.	Limite de temperatura
RO.	Limita de temperatură
AR.	حدود درجة الحرارة
HU.	Hőmérsékleti határértékek
FI.	Lämpötilarajoitus
NO.	Temperaturbegrensning
ID.	Batasan suhu

MD

CE

ZH. EO 滅菌
CS. Sterilní - Sterilizováno ethylenoxidem
KK. Стерильді - Этилен тотығымен стерильденген
PT. Estéril - Esterilizado por óxido de etileno
RO. Steril - Sterilizat cu oxid de etilenă
AR. معقم - معقم بأكسيد الإيثيلين
HU. Steril - Etilén-oxiddal sterilizált
FI. Steriili - steriloitu etyleenioksidilla
NO. Steril, etylenoksid-sterilisert
ID. Steril - Etilen oksida disterilkan

EN. Medical device
IT. Dispositivo medico
FR. Dispositif médical
DE. Medizinprodukt
ES. Producto sanitario
TR. Tıbbi cihaz
EL. Ιατρική συσκευή
NL. Medisch hulpmiddel
RU. Медицинское изделие
PL. Wyrób medyczny
SR. Medicinsko sredstvo
LT. Medicinos priemonė
SV. Medicinteknisk produkt
BG. Медицинско изделие
DA. Medicinsk udstyr
HR. Medicinski proizvod
ZH. 醫療器械
ZH. 医疗器械
CS. Zdravotnický prostředek
KK. Медициналық құрылғы
PT. Dispositivo médico
RO. Dispozitiv medical
AR. جهاز طبي
HU. Orvostechnikai eszköz
FI. Lääkinnällinen laite
NO. Medisinsk utstyr
ID. Perangkat medis

CS. Nepoužívejte, pokud je poškozený obal
KK. Орам зақымдалған болса, пайдаланбаңыз
PT. Não utilizar se a embalagem estiver danificada
RO. A nu se utiliza dacă ambalajul este deteriorat
AR. لا تستخدمه في حالة تلف الحزمة
HU. Ne használja, ha a csomagolás sérült
FI. Älä käytä, jos pakkaus on vaurioitunut
NO. Må ikke brukes hvis emballasjen er skadet
ID. Jangan digunakan jika kemasan rusak

EN. Double sterile barrier system
IT. Sistema a doppia barriera sterile
FR. Système à double barrière stérile
DE. Doppeltes Sterilbarrieresystem
ES. Sistema de doble barrera estéril
TR. Çift steril bariyer sistemi
EL. Σύστημα διπλού στείρου φραγμού
NL. Dubbel steriel barrièresysteem
RU. Двойная стерильная барьерная система
PL. System podwójnej bariery sterylnej
SR. Sistem sa dve sterilne barijere
LT. Dviguba sterilių barjerų sistema
SV. Dubbelt sterilit barriärsystem
BG. Двойна стерилна бариерна система
DA. Dobbelt steril barrieresystem
HR. Sustav dvostruke sterilne obloge
ZH. 双层无菌屏障系统
ZH. 双层无菌屏障系统
CS. Systém se dvěma sterilními bariérami
KK. Екі стерильді тосқауыл жүйесі
PT. Sistema de barreira estéril duplo
RO. Sistem de barieră sterilă dublă
AR. نظام عازل معقم مزدوج
HU. Kettős sterilgát-rendszer
FI. Kaksinkertainen steriili estojärjestelmä
NO. Dobbelt steril barrieresystem
ID. Double sterile barrier system

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■ ENGLISH ■

DESCRIPTION

Disposable device manufactured/assembled by RAND. Not made with natural rubber latex. Made of non-pyrogenic, non-toxic and biocompatible materials according to ISO 10993 standards and in relation with the intended use.

Specifically designed for use in oncology for intraperitoneal and intrapleural hyperthermic perfusion.

TECHNICAL FEATURES

Diameter: 8 mm

Perforated section: 230 mm

Total length: 600 mm

Use temperature: from 37 to 45°C

INTENDED USE

Surgically invasive medical device intended to infuse and drain clear solution also enriched with chemotherapeutic agents, in hyperthermic conditions, in thoracic or abdominal cavity. Disposable device also utilizable for surgical wound post-operative drainage.

SIDE EFFECTS AND CONTRAINDICATIONS

There is no evidence of side effects specifically related to the use of the product. The device is contraindicated for procedures outside of the intended use.

STORAGE CONDITIONS

Store the product in a temperature range 5°C – 35°C (41°F – 95°F).

INDICATIONS, EXPECTED CLINICAL BENEFITS, PATIENT TARGET GROUP, INTENDED USERS

Please refer to the Performer equipment user manual.

SAFETY INFORMATION

Information intended to attract the attention of the user to prevent potentially dangerous situations and to ensure correct and safe use of the device is indicated in the text in the following way:

WARNING

Indicates serious adverse reactions and potential safety hazards for the user and/or the patient that can occur both during proper use or misuse of the device, and also the limitations of use and the measures to be adopted in such cases.

CAUTION

Indicates any special care the user must exercise for the safe and effective use of the device.

For further information and/or in case of complaint contact RAND or the authorised local distributor.

WARNING

- The user must read and understand all user information prior to use. Failure to read and follow all instructions, or failure to observe all stated warnings, could cause serious injury or death to the patient.
- The device is intended for professional use. The procedure must be carefully and constantly monitored by qualified and trained personnel.
- Carry out a visual inspection and carefully check the devices before use. Transport and/or storage conditions other than those prescribed may have caused damage to the devices.
- Sterility is guaranteed only if the sterile packaging is not wet, opened, damaged or broken. Do not use the devices if sterility cannot be guaranteed.
- Check the expiry date on the label attached. Do not use the devices after the date shown.
- The devices must be used immediately after opening the sterile packaging.
- For single use only. Reuse of the devices brings with it the risk of cross-infection regardless of the cleaning or sterilization method used.
- Do not resterilise, do not undergo any further processing. Resterilization may not be effective and may compromise the integrity of the device thus resulting in unintended injury.
- Do not alter or manipulate the device in any way, in particular do not get out extra holes.
- Do not allow alcohol, anaesthetic fluids (i.e. isoflurane) or corrosive solvents (i.e. acetone) to come into contact with the device as they may jeopardize the structural integrity.
- The device is designed for a maximum working pressure of 500 mmHg; Rand is not responsible for improper setting of the parameters or accidental over-pressurization of the devices exceeding this safety limit.

CAUTION

- After use, dispose of the device according with the regulations in force in the country where the device is used.

SET UP (general information)

WARNING

The device must be handled aseptically.

1. Extract the devices from the external packaging.
2. If necessary perform a contra-incision of proper dimensions to allow insertion of the device.
3. Position the catheter through the contra – incision into selected seat.
4. Fix the device by means of a stitch or plaster.
5. Connect the catheter with the withdrawal or infusion line through the suitable asymmetrical connector placed on the line itself.

WARNING

The positioning and fastening procedures, the choice of the catheter types, as well as incision and anchorage method are under the direct responsibility of the surgical staff.

WARNING

The device is intended for body cavities use only, and not for vascular and/or orifice insertion.

CAUTION

1. During fastening carefully avoid stitching in order to prevent tearing or damaging of the device.
2. In case of anchorage with suture thread, make sure that drainage is not damaged by needles or sharp edges.
3. After positioning, check that the perforated section is completely included into cavity, moreover verify that the connection between line and catheter is perfectly sealed.
4. During the therapy check that devices placed on withdrawal lines are fairly triggered (the perforated portion must be completely immersed).

REMOVAL PROCEDURE

1. Remove the tape or the fastening stitch.
2. Disconnect the device from its respective line.
3. Remove the device by hand pulling slightly, avoid sudden movements. Do not use metal instruments.

CAUTION

If the device is left inside cavity for a long period of time, it might be difficult to remove it. Proceed cautiously.

WARNING

After removal ensure complete integrity of the device in order to exclude that any device fragments remain in the cavity.

WARNING

The devices are intended for use as drainage system up to a maximum of 24 hours (6 hours for Hyperthermic treatments with infusion/drainage of chemotherapeutic solutions). RAND does not assume any responsibility for an improper use, a wrong positioning/removal of the device and for the direct consequences for the patient and/or the operator in case the usage limits shown in this instruction for use are exceeded.

WARNING

The compatibility of use of the quoted devices with some commonly used medications - at the discretion of the physician - in intraperitoneal and/or intrapleural hyperthermal perfusion treatments has been appropriately verified by careful review of scientific literature and analysis of clinical test data relating to use of commercial medications based on the following active ingredients of consolidated use: CARBOPLATIN, CISPLATIN, MITOMYCIN C, ADRIAMYCIN/DOXORUBICIN, IRINOTECAN, MELPHALAN and OXALIPLATIN. Rand does not assume any responsibility for use of the devices in combination with medications or preparations not previously evaluated.

MEDICAL DEVICES FOR USE IN COMBINATION

RAND catheters are expressly designed for use in combination with infusion or withdrawal lines included in the following RAND products intended for Hyperthermic Intra-Peritoneal and/or Intrapleural Perfusion: HANG&GO 3 R9900141, HANG&GO 3 PLUS R9900140, HANG&GO HT R9900067 and HANG&GO HT BASIC R9900088.

REPORTING OF INCIDENT

Any serious incident that has occurred in relation to the device should be reported directly to Rand or to its distributor and the competent authority of the Member State in which the user is established.

WARNING

RAND does not assume any responsibility in case of use of alternative products other than those specified, unless they have been previously evaluated.

LIMITED WARRANTY

This limited warranty is in addition to any statutory rights of the purchaser pursuant to applicable law.

RAND warrants that all reasonable care has been taken in the manufacture of this medical device, as required by the nature of the device and the use for which the device is intended. RAND warrants that the medical device is capable of functioning as indicated in the current instruction for use when used in accordance with them by a qualified user and before any expiry date indicated on the packaging. However, RAND cannot guarantee that the user will use the device correctly, nor that the incorrect diagnosis or therapy and/or that the particular physical and biological characteristics of an individual patient, do not affect the performances and effectiveness of the device with damaging consequences for the patient, even though the specified instruction for use have been respected. RAND, whilst emphasizing the need to adhere strictly to the instruction for use and to adopt all the precautions necessary for correct use of the device, cannot assume any responsibility for any loss, damage, expense, incidents or consequences arising directly or indirectly from the improper use of this device. RAND undertakes to replace the medical device in the event that it is defective at the time of placing on the market or whilst being shipped by RAND up to the time of delivery to the final user unless such defect has been caused by mishandling by the purchaser. The above replaces all other warranties explicit or implicit, written or verbal, including warranties of merchantability and fitness for purpose. No person, including any representative, agent, dealer, distributor or intermediary of RAND or any other industrial or commercial organization is authorized to make any representation or warranty concerning this medical device excepted as expressly stated herein. RAND disclaims any warranty of merchantability and any warranty of fitness for purpose with regard to this product other than what is expressly stated herein. The purchaser undertakes to comply with the terms of this limited warranty and in particular agrees, in the event of a dispute or litigation with RAND, not to make claims based on alleged or proven changes or alterations made to this Limited Warranty by any representative, agent, dealer, distributor or other intermediary. The existing relations between the parties to the contract (also in the case that it is not drawn up in writing) to whom this Warranty, its interpretation and execution, nothing excluded and/or reserved, are regulated exclusively by the Italian law and jurisdiction. The court chosen is the Court of Modena (Italy).

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RanD S.p.A.